

## Product datasheet for **RC227044L3V**

### ADCK1 (NM\_001142545) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ADCK1 (NM_001142545) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADCK1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001142545
ORF Size:	1365 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227044).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001142545.1</a>
RefSeq ORF:	1368 bp
Locus ID:	57143
UniProt ID:	<a href="#">Q86TW2</a>
Cytogenetics:	14q24.3
Protein Families:	Druggable Genome, Protein Kinase, Secreted Protein
MW:	51.8 kDa
Gene Summary:	The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).[UniProtKB/Swiss-Prot Function]



[View online »](#)