

Product datasheet for **RC227053**

ALGIL2 (NM_001136152) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ALGIL2
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC227053 representing NM_001136152
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAGCTACTGCAGGCTGGGCTGTGACCGTCTACGACAAGCCGGCATCTTTCTTTAAAGAGGCACCTC
TGGACCTGCAGCACCGGCTCTTCATGAAGCTGGGCAGCAGCACTCTCCGTTCCAGGGCCCGCTCAGAACC
TGAGGACCCAGACAGAGCGGTCGGCCTTCACGGAGCGGGATTCTGGGAGCGGGCTGGTGACGCGTCTC
CAGGAGCGGCCAGCCCTGCTGGTCAGCAGCACAAAGCTGGACAGAGTTTGAACAAGTACTCTTGACGGAC
AGAACCTTCCTTCTCTCGTCTGTGTGATAACAGGCAAGGGCCTCTGAGGGAGTATTACAGCCGCTCAT
CCACCAGAAGCATTTCCAGCACATCCAGGTCTGCATCCCCTGGCTGGAGGGCCGAGGACTACCCCGCTT
CTAGGGTCGGTGGACCTGGATGTCTGTCTGGACACGTCTCCAGTGGCCTGGACCTGCCCATGAAGGTGG
TGGACATGTTCAAGGTGCTGTTTGCCTGCGTGTGCCGTGAAGTCAAGTGTTTACATGAGCTGGTGAAACA
TGAAGAAAACCGCCTGGTCTTTGAGGACTCAGAGGAAGTGGCAGCTCAGCTGCAGTATTTTGCAGATGCT
TTTCTCAAATTTCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227053 representing NM_001136152
Red=Cloning site Green=Tags(s)


MGATAGWAVTVYDKPASFFKEAPLDLQHLFMKLGSTHSPFRARSEPDPDTERSAFTERDSGSGLVTRL
HERPALLVSSTSWTEFEQLTLDGQNLPSLVCVITGKGPLREYYSRLIHQKHFIQVCIPIWLEGRGLPPL
LGSVDLDVCLDTSSSLDLPKVVDMFRCLPACAVNFKCLHELKHEENRLFEDSEELAAQLQYFADA
FLKLS

SGPTRRRLE**QKLISEEDLA**AND**I**LDYKDDDDKV



Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001136152

ORF Size: 645 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:	NM_001136152.1, NP_001129624.1
RefSeq ORF:	648 bp
Locus ID:	644974
UniProt ID:	C9J202
Cytogenetics:	3q22.1
Protein Families:	Transmembrane
MW:	24 kDa
Gene Summary:	Putative glycosyltransferase.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC227053



Western blot validation of overexpression lysate (Cat# [LY427825]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227053 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).