

Product datasheet for **RC227067**

BHLHB9 (NM_001142530) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | BHLHB9 (NM_001142530) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | BHLHB9 |
| Synonyms: | GASP3; p60TRP |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>RC227067 representing NM_001142530
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCTGGGACTAAGAATAAGACAAGAGCCAGGCCAAAAGCTGAAAAAAGGCTGCTATACAAGCTAAAG
 CTGGAGCAGAGAGGGAGGCTACTGGTGTGTTAGGCCTGTAGCCAAGACCAGGGCCAAAGCAAAGCCAA
 GACAGGGTCTAAGACAGATGCAGTAGCAGAGATGAAGGCAGTGTCTAAGAACAAGGTTGTTGCTGAGACG
 AAGGAAGGAGCTCTGTCAGAGCCTAAGACTCTGGGCAAAGCCATGGGAGATTTCACTCCCAAGGCTGGGA
 ATGAGTCCACCAGCTCCACATGTAAAAATGAGGCTGGTACTGATGCCTGGTTCTGGGCTGGGGAAGAGGC
 CACTATCAATTCCTGGTCTGGAATGGAGAAGAGGCTGGTAATAGTTTCAGCACTAAGAATGATAAACCT
 GAAATTTGGTGGCCAGGCTGTGCTGAGGAGTTGGAACCTGCGGCTGGGGCCGATTGCAAACCTAGGTCAG
 GGGCTGAGGAGGAGGAGGAAGAGAATGTTATTGGGAAGTGGTTTGGGAAGGAGATGATACTAGTTTTGA
 CCCTAATCCTAAACCTGTGAGCAGGATAGTTAAGCCTCAGCCTGTGTATGAAATTAATGAAAAAATAGG
 CCCAAGGACTGGTCTGAGGTAATCTGGCCCAATGCCCTGCTGTAACCTCCAGCTGTGTTAGGATTTA
 GATCACAGGCACCATCTGAGGCAAGCCCTCCTTCATATATTGTTCTGGCCTCCGCTGAAGAAAATGCCTG
 TTCTTTGCCTGTGGCAACAGCTTGCCGCCCTTCTAGGAACACTCGCTCATGCTCACAGCCTATCCCTGAG
 TGTCTGTTTTGATTCTGACCCCTGCATCCAGACCATAGATGAGATTAGACGTCAAATCAGGATCAGGGAGG
 TAAATGGGATTAAGCCATTTGCTGTCTTGCAAAATGGAATGCTATATGGATTCTGAGGAATTTGAAAA
 ACTTGTAGCTTACTTAAGTCACTACTGATCCTCTTATTCATAAAAATAGCACGGATTGCAATGGGTGTC
 CATAATGTTCACCATTTGCCAAAGAGTTTATTAACGAAGTAGGTGTAGTGACACTTATTGAAAGCTTGC
 TCAGTTTTCTTCCCCTGAAATGAGAAAAAGACTGTAATTAAGTCTGAATCCTCCTTCTGGGGATGAAAG
 ACAACGCAAAAATTGAATTACATGTTAAGCATATGTGTAAGAAAACCATGTCATTTCTTTGAACTACCG
 GGACAGCAATCTGGATTAAGATACTAGGACAAGTACTGACTACTGATTTTGTCCATCACTACATTGTTGCCA
 ATTACTTTTCAGAGCTTTTCCATTTGCTGTCTCAGGAAATGCAAAAACAGAAATCTTGTTTTGAAGT
 ACTTTTAAATATGCTGAAAATCCAAGTGCAGCCAGAGACATGATCAATATGAAGGCATTGGCAGCATTAA
 AAAGTCACTTTAACCAGAAAGAGGCAAAAGCCAACTTGTAGTGGTGTGGCCATTTTATTAACATAA
 AGGAGCATATCAGAAAAGGCTCAATTGTAGTTGTTGATCACTTGAGTTATAATACACTCATGGCCATTTT
 CAGGGAAGTTAAGAGATTATTGAAACAATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC227067 representing NM_001142530
 Red=Cloning site Green=Tags(s)

MAGTKNKTRAQAKTEKKAIIQAKAGAEREATGVVRPVAKTRAKAKAKTGSKTDAVAEMKAVSKNKVVAET
 KEGALSEPRTLKAMGDFTPKAGNESTSSCKNEAGTDAWFWAGEEATINSWFWNGEEAGNSFSTKNDKP
 EIGAQVCAEELEPAAGADCKPRSGAEEEEENVIGNWFWEGDDTSFDPNPKPVSRIKVPQPVYEINEKNR
 PKDWSEVTIWPNAVTPAVLGFERSQAPSEASPPSYIVLASAENACSLPVATACRPSRNRTRSCSQPIPE
 CRFSDSPCIQTIDEIRRQIRIREVNGIKPFACPKMECYMDEEFELVSLKSTTDPLIHKIARIAMGV
 HNVHPFAQEFINEVGVVTLIESLLSFPSPEMRKTIVITLNPPSGDERQRKIELHVHKMKETMSFPLNSP
 GQQSGLKILGQLTTDFVHHYIVANYFSELFHLLSSGNCKTRNLVLKLLNMSNP TAARDMINMKALAAL
 KLIFNQEKANLVSGVAIFINIKHEIRKGSIVVVDHLSYNTLMAIFREVKEI IETM

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001142530

ORF Size: 1641 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142530.1](#), [NP_001136002.1](#)

RefSeq ORF: 1644 bp

Locus ID: 80823


UniProt ID: [Q6PI77](#)

Cytogenetics: Xq22.1

MW: 60.1 kDa

Gene Summary: This gene is a member of a gene family which encodes proteins with a basic helix-loop-helix domain. Other members of this gene family encode proteins which function as transcription factors, either enhancing or inhibiting transcription depending on the activity of other DNA binding proteins. The coding region of this gene is located entirely within the terminal exon. The encoded protein may be involved in the survival of neurons (PMID: 15034937). Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2011]

Product images:



Circular map for RC227067