

Product datasheet for **RC227072L3V**

SCRN1 (NM_001145513) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SCRN1 (NM_001145513) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCRN1
Synonyms:	SES1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145513
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227072).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145513.1 , NP_001138985.1
RefSeq Size:	5286 bp
RefSeq ORF:	1245 bp
Locus ID:	9805
UniProt ID:	Q12765
Cytogenetics:	7p14.3
MW:	46.4 kDa



[View online »](#)

Gene Summary:

This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]