

Product datasheet for **RC227075**

AGAP5 (NM_001144000) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	AGAP5
Synonyms:	CTGLF2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC227075 representing NM_001144000
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGGAACATACTGACCTGTTGTGTGCACCCTAGCGTCAGCCTCGAGTTTGACCAGCAACAGGGGTCGG
 TGTGTCCCTCTGAATCTGAGATCTATGAGGCAGGAGCTGGGGACAGGATGGCAGGAGCGCCATGGCTGC
 TGCTGTGCAGCCTGCTGAGGTGACCGTTGAAGTTGGTGAGGACCTCCACATGCACCACATTCGTGACCAG
 GAGATGCCTGAAGCTTTGGAGTTTACCCCTTCTGCCAATCCAGAGGCAAGCACAATATTCAGAGGAACTC
 TCAAAACAGATGCTTTGGAGTTTAAACCCTTCTGCCAATCCAGAGGCAAGCACAATATTCAGAGGAACTC
 TCAAAACAGATGTTGTAGAAATAGAAGAAGCAACTGTACAAACCATGTATCTACTGAGCGTTTCAGTCAA
 CAATACAGCTCGTGTTCGACAATATTCCTTGATGACAGCACAGCCAGCCAGCATTATCTTACAATGACAA
 TAATATCTGTGACCTTGAGATACCTCATCATATCACACAAAGAGATGCAGATAGATCTTTGAGCATACC
 TGATGAACAGTTACACTATTTGCGGTTTCCACCGTGACATTACGAAGAACAGAAATGGAGGTGGGAGT
 TTAAATAACTATTCTCTCCATTCCATCGACTCCCAGCACCAGCCAGGAGGACCCTCAGTTCAGTGTTT
 CTCCCCTGCAACACACCCACCCCGTTTGCAAGCGGTCCATGCGCTGGTCCAACCTGTTTACATCTGA
 GAAAGGGAGTCACCCAGACAAAGAGAGGAAAGCCCGGAGAATCATGCTGACACCATCGGGAGCGGCAGA
 GCCATCCCCATTAAACAGGGCATGCTCTTAAAGCGAAGTGGGAAATGGCTGAAGACATGAAAAAGAAAT
 ACGTCACCGTGTGTTCCAATGGCGTGCTCACCTATTATTCAAGCTTAGGTGATTATATGAAGAAATTTCA
 TAAAAAGAGATTGACCTTCGACATCTACCATCAAAGTCCCAGGAAAGTGGCCATCCCTAGCCACATCG
 GCCTGTGCACCCATCTCCAGCTCTAAAAGCAATGGCCTATCCAAGGACATGGACACCGGGTGGTGACT
 CCATATGCTTCAGCCCCAGTATCTCCAGCACCACCAGCCCCAAGCTCAACCCGCCCCCTCTCTCATGC
 CAATAAAAAGAAACACCTAAAGAAAGAAAGCACCACCAACACTTTATGATTGTGTCTGCCACTGGCCAAACA
 TGGCACTTTGAAGCCACGACGTATGAGGAGCGGGATGCCTGGGTCCAAGCCATCCAGAGCCAGATCCTGG
 CCAGCCTGCAGTCATGCGAGAGCAGTAAAAGCAAGTCCCAGCTGACCAGCCAGAGCGAGGCCATGGCCCT
 GCAGTCGATCCAAACATGCGTGGGAACGCCCCTGTGTGGACTATGAGACCCAGAATCCTAAGTGGGCC
 AGTTTGAAGTGGGAGTCTCATGTGTATTGAATGCTCAGGAATCCACCGCAGTCTTGGCACCCGCTTT
 CCCGTGTGCGATCTCTGGAGCTGGATGACTGGCCAGTTGAGCTCAGGAAGGTTATGTCATCTATTGGCAA
 TGACCTAGCCAACAGCATCTGGGAAGGGAGCAGCCAGGGGCAGACAAAACCTCAGTAAAGTCCACGAGG
 GAAGAGAAGGAACGGTGGATCCGTTCCAAATATGAGGAGAAGCTCTTTCTGGCCCCACTACCCTGCACTG
 AGCTGTCCCTGGGCCAGCACCTGCTGCGGGCCACCGCTGATGAGGACCTGCAGACAGCCATCCTGCTGCT
 GGCACATGGCTCCCGTGAGGAGGTGAACGAGACCTGTGGGGAGGGAGACGGCTGCACGGCGCTCCATCTG
 GCCTGCCGCAAGGGGAATGTGGTCTGGCAGACCTCTGATCTGGTACGGGGTGGACGTATGGCCCCGAG
 ATGCCACGGGAACACAGCGCTGACCTACGCCCGCAGGCCCTCCAGCCAGGAGTGCATCAACGTGCTTCT
 GCAGTACGGCTGCCCGACGAGTGCGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227075 representing NM_001144000
 Red=Cloning site Green=Tags(s)


MGNILTCVHPSVSLFDDQQQSVCPSSESEIYEAGAGDRMAGAPMAAAVQPAEVTVEVGEDLHMHIRDQ
 EMPEALEFSLSANPEASTIFQRNSQTDALFNPSANPEASTIFQRNSQTDVVEIRRSNCTNHVSTERFSQ
 QYSSCSTIFLDDSTASQHYLTMTIISVTLEIPHHITQRDADRSLIPDEQLHSFAVSTVHITKNRNGGGS
 LNNYSSSIPSTPSTSQEDPQFSVPPTANTPTPVCKRSMRWSNLFTSEKGSHPDKERKAPENHADTIGSGR
 AIPKQGMILLKRSGLKLTWKKKYVTLCSNGVLTYYSSLGDYMKNIHKKEIDLRTSTIKVPKWPPLATS
 ACAPISSSKSNGLSKMDTGLGDSICFSPISSTTSPKLNPPPSPHANKKKHLKKKSTNNFMIVSATGQT
 WHFEATTYEERDAWVQAIQSQILASLQSCSSKSKSLTSQSEAMALQSIQNMGRNAHCVDYETQNPKW
 SLNLGLVLMCIECSGIHRSGLTRLRVSRLSLDDWPVELRKVMSSIGNDLANSIWEQSSQGQTKPSVKSTR
 EEKERWIRSKYEELFLAPLPCTELSLGQHLLRATADEDLQTAILLAHGSREEVNETCGEGDGCTALHL
 ACRKGNVLAQLLIWYGVDMARDAHGNTALTARQASSQECINVLLQYGCPDECV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8033_a05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001144000

ORF Size: 2058 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22µm filter is required.

RefSeq:	NM_001144000.2
RefSeq ORF:	2061 bp
Locus ID:	729092
UniProt ID:	A6NIR3
Cytogenetics:	10q22.2
MW:	75.3 kDa
Gene Summary:	Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC227075



Western blot validation of overexpression lysate (Cat# [LY428460]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227075 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).