

Product datasheet for **RC227075L3V**

AGAP5 (NM_001144000) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	AGAP5
Synonyms:	CTGLF2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001144000
ORF Size:	2058 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC227075).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001144000.1
RefSeq ORF:	2061 bp
Locus ID:	729092
UniProt ID:	A6NIR3
Cytogenetics:	10q22.2
MW:	75.3 kDa



Gene Summary: Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]