

Product datasheet for RC227090

FAM185A (NM_001145268) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FAM185A
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC227090 representing NM_001145268
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTTGCCCCCTGCTCAGGTTGGGAGCTTGGCTGCTTCCGTCTGTCTCCGTCAGGTCCGACTGTGGG
CTGGCGCTGGGCGCTGGGCTTGCTGGGCTTGCCAAGCCAGGCCGTACAGCTCAGGTGGGAGCGAGCGCTG
GCCCCGATCGGAGACTGAGGTCCCTCCGCTGGCCGCGCGCCGAACCTGAAGGAGTGGACACTGCAG
GTGAGCCCCGTTTGGTGGCTGCGGGCGCGGCTCCCGTGCCACCTGGCCGTGAGGCCCTGGACCCCTCA
CCTACCCGGATGGCGACCGCGTGTGGTGGCGGTGCGGCGTGGAGGGCGCGTGGGGGCTGGACGG
CCTGCAGGTGAAGTACGACGAGGATCTGGAGGAGATGGCCATTGTGTCTGATACTATCCACCCCGAGCG
TCCGTGGAGGTGAACGCGCCCTGAAGTTTGGTTTAGATATCAAGTCATCAGGGTCTGGCTGTGTAAG
TTCAAAGTATTGAGGTGATAATTGCAAAATTGAAACAGAGCATGGGACTAGTATCTTGCACTGTGTTAA
GGGTCAAAATTGCATGTTCAAACAAAGGAGGCAAAGTGATCTGTCTGGGAACAGTTTATGGAATATA
GATATTCATGCATCAGATAAAAGTGTGTGACCATAGATAAACTGCAGGGAAGTTCTGTTACTGTATCTA
CCGAAGATGGTTTGTGAAAGCCAAGTATCTTTATACAGAATCATATTTCTGTCTTCTGCTGCTGGGA
TATTACATTAGGAAGTGTTCATGTAATATAACATTACAAAGCAAGATGGGTAACATCACAGTAGATTCTG
TCTTCTGGATGTCTAAAGCCTCAACTAATCAGGGTGCCATAGATGTTTATGTCAGCCAACCTGGGAAAG
TGGAATTGAAATCCCATAAAGGTTCTATTATTGTCAAGGTCCCCTCTTCACTTCAAGCTCATTTACAGTT
ATCAGGGAAGAGGTTGATGTGAAGTCAAGTCCATGTTTCAAGAAATGGCTGAAGTTCTGTAAGATGAT
GTTGTAACAGTGACTGGACTCATGAATCAAGCAAGCAACGTGAAAAATGGATTAAGGCTGATGCCCAA
AAGGCACTGTAAGTTTGAAGGCAAAGTTGGTTTCAGTCCCTGAACTGCAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATACGACGATAAGGTTTAA



Protein Sequence:

>RC227090 representing NM_001145268

Red=Cloning site Green=Tags(s)

MLAPCSGWELGCFRLCLRQVRLWAGAGRWACQARPYSSGGSERWPGSETVPPPGPARRTLKEWTLQ
VSPFGRRLRRLPCHLAVRPLDPLTPDGDRLVAVCGVEGGVRGLDGLQVYDEDELEMAIVSDTIHPQA
SVEVNAPLKFGLDIKSSSGSCVKVQSIEGDNCKIETEHGTSILQSVKGQKLHVQTKGGKVICLGTVYGN
DIHASDKSAVTIDKLQSSSVTVSTEDGLLKAKYLYTESSFLSSAAGDITLGSVHGNITLQSKMGNITVDS
SSGCLKASTNQGAIDVYVSQLGKVELKSHKGSIIKVPSSLQAHQLSGKEVDVNSEVHVQEMAEVRKDD
VYTVTGLMNQASKREKWKADAPKGTVSFRRQSWFQSLKLQD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8033_a01.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001145268

ORF Size:

1176 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145268.1](#), [NP_001138740.1](#)

RefSeq ORF: 1179 bp


Locus ID: 222234

UniProt ID: [Q8N0U4](#)


Cytogenetics: 7q22.1

MW: 42.1 kDa

Product images:



Circular map for RC227090



Western blot validation of overexpression lysate (Cat# [LY428775]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227090 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).