

Product datasheet for **RC227108L3V**

KREMEN1 (NM_032045) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	KREMEN1
Synonyms:	ECTD13; KREMEN; KRM1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032045
ORF Size:	1476 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC227108).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_032045.3
RefSeq ORF:	1479 bp
Locus ID:	83999
UniProt ID:	Q96MU8
Cytogenetics:	22q12.1
Domains:	KR, CUB, WSC



Protein Families: Druggable Genome, Transmembrane

MW: 51.9 kDa

Gene Summary: This gene encodes a high-affinity dickkopf homolog 1 (DKK1) transmembrane receptor that functionally cooperates with DKK1 to block wingless (WNT)/beta-catenin signaling. The encoded protein is a component of a membrane complex that modulates canonical WNT signaling through lipoprotein receptor-related protein 6 (LRP6). It contains extracellular kringle, WSC, and CUB domains. Alternatively spliced transcript variants encoding distinct isoforms have been observed for this gene. [provided by RefSeq, Jul 2008]