

Product datasheet for RC227111

BRP44 (MPC2) (NM_001143674) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: Myc-DDK
Symbol: BRP44
Synonyms: BRP44; SLC54A2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC227111 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCGGCCGCGGTGCCCGAGGCCTGCGGGCCACCTACCACCGGCTCCTCGATAAAGTGGAGCTGATGC
TGCCCCGAGAAATTGAGGCCGTTGTACAACCATCCAGCAGGTCCAGAACAGTTTCTTCTGGGCTCCAAT
TATGAAATGGGGGTTGGTGTGTGCTGGATTGGCTGATATGGCCAGACCTGCAGAAAACTTAGCACAGCT
CAATCTGCTGTTTGTATGGCTACAGGGTTATTTGGTCAAGATACTCACTTGTAAATTATCCAAAAAATT
GGAGTCTGTTTGTCTGTTAATTTCTTTGTGGGGCAGCAGGAGCCTCTCAGCTTTTCTGTAATTGGAGATA
TAACCAAGAACTAAAAGCTAAAGCACACAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227111 protein sequence
Red=Cloning site Green=Tags(s)

MSAAGARGLRATYHRLLDKVELMLPEKLRPLYNHPAGPRTVFFWAPIMKWGLVCAGLADMARPAEKLSTA
QSAVLMATGFIWSRYSLVIIPKNWSLFAVNFFVGAAGASQLFRIWRYNQELKAKAHK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6457_d11.zip

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001143674

ORF Size: 381 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001143674.4](#)

RefSeq Size:	2220 bp
RefSeq ORF:	384 bp
Locus ID:	25874
UniProt ID:	O95563
Cytogenetics:	1q24.2
MW:	14.3 kDa
Gene Summary:	Mediates the uptake of pyruvate into mitochondria.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC227111



Western blot validation of overexpression lysate (Cat# [LY428311]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227111 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).