

Product datasheet for **RC227113**

ACOXL (NM_001142808) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ACOXL
Synonyms:	acyl-Coenzyme A oxidase-like; FLJ11042
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC227113 representing NM_001142808
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAGAGCTTTGACAGTTCAGAGAGTGAAGTTTGCCATGGACCTGCCTCTGTTAAACGTGCAGGTCAGG
 ATCTGGCAGAGAAAACAAAGATTTTGTACGCCGAAGCCTTGTATAGGAGAAGTCCTCTCCATGGCGGA
 CATGGCCACAGGAGTGAAGTGCAGGAATAATTTATTGGCTATTTGGTGGTCTATCAGGAATCTCGGAAGC
 CCTGAACATGTTACTAAGTGGTTTCAGCCACTCCAGGAGCAGAAATACACTGGGATGTTTGAATGACCG
 AGAGGGGCCATGGGAGCAACGCGAGAGGGATCCAGACCGAAGCCACCTTTGACCTCTCTGCCAGGAGTT
 TGTAAATTGACACGCCGTGTGAAAATGCGGAGAAGATGTATATTGGAAATGCCATGTACGGGAATTTATGCA
 GCTGTCTTTGCCAGCTCATCATAGATGGAAGATCTCAAGGGCCCCACTGTTTCATCGTTCCTGTCCGGG
 ATGAAAACGGAAGCTTGATCCAGGAGTCACAGCTATTGATATGATGTACAAGGAGGGTCTGCATGGTGT
 GGACAAATGGGATATTAATTTTGACAAGTTTCGGATACCCAGGAGAACCTGCTGGATAAGTTTGGTTCC
 GTGGCTCCAGATGGACAGTACCATTGCGCTATTAGGAACAAGAGTGCAAGATTCAATGCCATGCTGGCAG
 CACTGACCCCTTCGAGATTAGCTGTGGCTTTCCAAGCTATGGGTGCCATGAAGCTTGGGTTGACGATAGC
 CATTGCTATAGCCACAGCCGAGGAGGCTTTGGGCCAAAACCAAGGAAGAGGTGAAGATCATTGAGCAC
 CAAACACAGACCTGCGGCTGATGCCCCACCTGGCCACAGCCTTGCCCTGACCTTCGTGACGAGGTATG
 CTGGGGCCCTCCTGGATGAGGATGTCTTCCAGGAAAGGAGCTGGTCAACAGTCGCTCGCTGCAGGCTCT
 GGTGGCGGGCTGAAGGCCTACAGCACCTGGGAGAACATCCGCTGCCTGCAGGACTGCCGCGAGTCACT
 GGAGGCATGGTTGTGGGCGGGAAGTGTGGCCAATACACAAACAGTATGAAGAAAACCACTCTTTG
 GCCTGCTCCAAACTGGGCTGAATCTGTGGGGACAAGCTGAGAACAGTTTCTTGGCATTAAACATGGA
 CACAGTTGATGATCTCGCTTTCTGTTGAAAGCAGTGAATTTCTGTAAAGGGTTCTTCAGCGGGGTTTG
 GTGGCCAGAATTTATTATAAGGTAAAGACCAAGAAGGAGGATTTTTCCATGCCTGGAACCTCGTGTCTGC
 ACCACGTGGCTTCTGTCCCTGGCACACACTACCGAGTTACCTTAGAGCAGTTCTCCCTAGCAGTGAA
 GAGCTGTCTGACCAAGAGGACAGACTTTGTTAATGAAGTTTGTCTGTTGTATGGAACCAAGCTGGTG
 TTTCAGGAGCGGGCTGGTATTTAGAACATAAACTTGACTCCCATGGCCAGCACGAGGATCAGGAATC
 AGGAGAGATGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227113 representing NM_001142808
 Red=Cloning site Green=Tags(s)


MRALTVQRVKFAFDLPLKRAQDLAEKTKNFVSRSLVIGEVLSMADMATGVKCGIIYWLFGGAIRNLGS
 PEHVTKWFQPLQEYKTMFAMTERGHGSNARGIQTEATFDLSAQEFVIDTPCENAEKMYIGNAMYGNYA
 AVFAQLIIDGRSQPHCFIVPVRDENGSLYPGVTAIDMMYKEGLHGVDNGILIFDKVRIPRENLLDKFGS
 VAPDQYHSPIRNKSARFNAMLAALTPSRLAVAFQAMGAMKLGTLIAIRYSHSRQFGPKTKEEVKIEH
 QTQTLRLMPHLATALALTFVSRYAGALLDEDVFQKELVNSRLQALVAGLKAYSTWENIRCLQDCRECT
 GGMVVGRELLAQYTKQYEEKPLFGLLQNWAEVSGDKLRTSFLAFNMDTVDDLAFLLKAVKFRERVLQRGL
 VARIYYKVKTKKEDFFHAWNSCLHHVASLSLAHTRVTLQFSLAVKSCPDQEDQTLMLKFCLLYGTCLV
 FQERAWYLEHKYLTMASTRIRNQCRC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1409_e11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142808

ORF Size: 1551 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001142808.1](#), [NP_001136280.1](#)


RefSeq ORF: 1553 bp

Locus ID: 55289
Cytogenetics: 2q13
MW: 58.2 kDa

Product images:



Circular map for RC227113



Western blot validation of overexpression lysate (Cat# [LY428287]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227113 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).