

Product datasheet for **RC227121L4V**

iASPP (PPP1R13L) (NM_001142502) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	iASPP (PPP1R13L) (NM_001142502) Human Tagged ORF Clone Lentiviral Particle
Symbol:	iASPP
Synonyms:	IASPP; NKIP1; RAI; RAI4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142502
ORF Size:	2484 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227121).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001142502.1 , NP_001135974.1
RefSeq Size:	3122 bp
RefSeq ORF:	2487 bp
Locus ID:	10848
UniProt ID:	Q8WUF5
Cytogenetics:	19q13.32
Protein Families:	Transcription Factors
MW:	89.1 kDa



[View online »](#)

Gene Summary:

IASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53; MIM 191170), whereas ASPP1 (MIM 606455) and ASPP2 (MIM 602143) are activators of p53.[supplied by OMIM, Mar 2008]