

Product datasheet for **RC227169L3V**

SPG20 (SPART) (NM_001142296) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SPG20 (SPART) (NM_001142296) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SPART |
| Synonyms: | SPG20; TAHCCP1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001142296 |
| ORF Size: | 1998 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC227169). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001142296.1 , NP_001135768.1 |
| RefSeq Size: | 5014 bp |
| RefSeq ORF: | 2001 bp |
| Locus ID: | 23111 |
| UniProt ID: | Q8N0X7 |
| Cytogenetics: | 13q13.3 |
| MW: | 72.9 kDa |



[View online »](#)

Gene Summary:

This gene encodes a protein containing a MIT (Microtubule Interacting and Trafficking molecule) domain, and is implicated in regulating endosomal trafficking and mitochondria function. The protein localizes to mitochondria and partially co-localizes with microtubules. Stimulation with epidermal growth factor (EGF) results in protein translocation to the plasma membrane, and the protein functions in the degradation and intracellular trafficking of EGF receptor. Multiple alternatively spliced variants, encoding the same protein, have been identified. Mutations associated with this gene cause autosomal recessive spastic paraplegia 20 (Troyer syndrome). [provided by RefSeq, Nov 2008]