

Product datasheet for **RC227183L3V**

Myotilin (MYOT) (NM_001135940) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Myotilin (MYOT) (NM_001135940) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Myotilin
Synonyms:	LGMD1; LGMD1A; MFM3; TTID; TTOD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135940
ORF Size:	1497 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227183).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001135940.1 , NP_001129412.1
RefSeq Size:	1808 bp
RefSeq ORF:	945 bp
Locus ID:	9499
UniProt ID:	Q9UBF9
Cytogenetics:	5q31.2
MW:	55.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a cytoskeletal protein which plays a significant role in the stability of thin filaments during muscle contraction. This protein binds F-actin, crosslinks actin filaments, and prevents latrunculin A-induced filament disassembly. Mutations in this gene have been associated with limb-girdle muscular dystrophy and myofibrillar myopathies. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.[provided by RefSeq, Oct 2008]