

Product datasheet for **RC227229L4V**

SLC25A44 (NM_001135672) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC25A44 (NM_001135672) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC25A44
Synonyms:	FLJ90431; KIAA0446; OTTHUMP00000016580; RP11-54H19.3; solute carrier family 25, member 44
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001135672
ORF Size:	615 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227229).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001135672.1 , NP_001129144.1
RefSeq ORF:	617 bp
Locus ID:	9673
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Transmembrane
MW:	23.3 kDa
Gene Summary:	SLC25A44 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).[supplied by OMIM, Mar 2008]



[View online »](#)