

Product datasheet for **RC227234L4V**

DIS3L (NM_001143688) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DIS3L (NM_001143688) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DIS3L
Synonyms:	DIS3L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001143688
ORF Size:	3162 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227234).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001143688.1
RefSeq ORF:	3165 bp
Locus ID:	115752
UniProt ID:	Q8TF46
Cytogenetics:	15q22.31
MW:	120.6 kDa


[View online »](#)

Gene Summary:

The cytoplasmic RNA exosome complex degrades unstable mRNAs and is involved in the regular turnover of other mRNAs. The protein encoded by this gene contains 3'-5' exoribonuclease activity and is a catalytic component of this complex. [provided by RefSeq, May 2016]