

## Product datasheet for **RC227243L4V**

### Invadolysin (LMLN) (NM\_001136049) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Invadolysin (LMLN) (NM_001136049) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Invadolysin
Synonyms:	GP63; INV; IX14; LMNL1; MSP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001136049
ORF Size:	2076 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227243).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001136049.2</a> , <a href="#">NP_001129521.2</a>
RefSeq ORF:	2055 bp
Locus ID:	89782
UniProt ID:	<a href="#">Q96KR4</a>
Cytogenetics:	3q29
Protein Families:	Druggable Genome, Protease, Transmembrane
MW:	77.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a zinc-metallopeptidase. The encoded protein may play a role in cell migration and invasion. Studies of a similar protein in *Drosophila* indicate a potential role in mitotic progression. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]