

Product datasheet for **RC227249L3V**

TYW1B (NM_001145440) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TYW1B (NM_001145440) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TYW1B
Synonyms:	LINC00069; NCRNA00069; RSAFD2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145440
ORF Size:	2004 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227249).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145440.1 , NP_001138912.1
RefSeq ORF:	2007 bp
Locus ID:	441250
UniProt ID:	Q6NUM6
Cytogenetics:	7q11.22-q11.23
MW:	76.7 kDa


[View online »](#)

Gene Summary:

Wybutosine is a hypermodified guanosine found in phenylalanine tRNA. Wybutosine functions to stabilize codon-anticodon interactions during ribosome decoding and therefore supports the maintenance of the reading frame. In yeast, the homolog of this gene is essential for the synthesis of wybutosine. The human genome contains two closely related genes that putatively function in wybutosine synthesis. The open reading frame of this locus is disrupted in some individuals. Thus, this locus appears to be an evolving pseudogene, but may still be functional in some members of the population. [provided by RefSeq, Apr 2014]