

## Product datasheet for **RC227269A1V**

### Human FXYD3 (NM\_001136012) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human FXYD3 (NM\_001136012) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** FXYD3  
**Synonyms:** MAT8; PLML  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC227269 representing NM\_001136012  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCAGAAGGTGACCCTGGGCCTGCTTGTGTTCTGGCAGGCTTTCCTGTCTGGACCCAATGACCTAG  
 AAGATAAAACAGTCCTTTCTACTATGACTGGCACAGCCTCCAGGTTGGCGGGCTCATCTGCGCTGGGGT  
 TCTGTGCGCCATGGGCATCATCATCGTCATGAGTGAGTGGAGGAGCTCGGGGAGCAGGCGGGCCGGGGC  
 TGGGGCTCCCCCTCCCCTGACCACTCAGCTCTCCCCAACAGGTGCAAAATGCAAAATGCAAGTTTGCCAGA  
 AGTCCGGTCACCATCCAGGGGAGACTCCACCTCATCACCCAGGCTCAGCCCAAAGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC227269 representing NM\_001136012  
 Red=Cloning site Green=Tags(s)

MQKVTLGLLVFLAGFPVLDANDLEDKNSPFYDWHSLQVGLICAGVLCAMGIIIVMSEWRSSGEQAGR  
 WGSPLTTQLSPTGAKCKCKFGQKSGHHPGETPPLITPGSAQS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001136012  
**ORF Size:** 339 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_001136012.1</u>
<b>RefSeq ORF:</b>	342 bp
<b>Locus ID:</b>	5349
<b>UniProt ID:</b>	<u>Q14802</u>
<b>Cytogenetics:</b>	19q13.12
<b>MW:</b>	11.99 kDa