

Product datasheet for **RC227271**

TSEN2 (NM_001145394) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	TSEN2
Synonyms:	PCH2B; SEN2; SEN2L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC227271 representing NM_001145394
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGCAGAAGCAGTTTTCCATGCCCAAAGAGGAAAAGAAGAGTGTATGAGACTTACGAGTCTCCATTGC
CAATCCCTTTTGGTCAGGACCATGGTCTCTGAAAGAATCAAGATATTCGGTCTGAAATGATTAAACA
CAATGTGATTGTGAGGAATGCGGAGGACATTGAGCAGCTCTATGGGAAAGGTTATTTGGAAAAGGTATT
CTTTCAAGAAGCCGTCCAAGCTTCAACAATTCAGATCCTAACTGGTTGCTAAATGGAAAGATATGAAGA
CAAACATGCCTATCATCACATCAAAGAGGTATCAGCATAGTGTGAGTGGGCAGCAGAGCTGATGCGTAG
ACAGGGGCAGGATGAGAGTACAGTGCGCAGAATCCTCAAGGATTACACGAAACCGCTTGAGCATCCTCCT
GTGAAAAGGAATGAAGAGGCTCAAGTGCATGACAAGCTTAACCTGGAATGGTTTCCAACATGGAAGGCA
CAGCAGGGGAGAGAGACCTTCTGTGCGTGGCCTTCATCATGAAGACGGCAGCCAGCACATCGGCCCTCT
GCATCCTGGGGACAGAGGGCCTGACCATGAGTACGTGCTGGTCGAGGAAGCGGAGTGTGCCATGAGCGAG
AGGGAGGCTGCCCCAAATGAGGAATTGGTGCAAAGAAACAGGTTAATATGCAGAAGAAATCCATATAGGA
TCTTTGAGTATTTGCAACTCAGCCTAGAAGAGGCCTTTTCTTGGTCTATGCTCTGGGATGTTTAAGTAT
TTACTATGAGAAGGAGCCTTTAACGATAGTGAAGCTCTGGAAGCTTTCACTGTAGTTCAGCCACGTTT
AGAACCACCTACATGGCCTACCATTAATCTTTCGAAGCAAGGGCTGGGTGCCCAAAGTGGGACTCAAGTACG
GGACAGATTTACTGCTATATCGGAAAGGCCCTCATTATACCATGCAAGTTATTCTGTCTATTATCGAGCT
AGTTGATGACCATTTTGAAGGCTCTCTCCGAGGCCTCTCAGTTGGAAGTCCCTGGCTGCCTTGAGCAGA
GTTTCCGTTAATGTCTTAAGGAATTATGCTGTGCTATTTGATTAAACCCTCTACTATGACTGACAAGG
AAATGGAGTCACCAGAATGTATGAAAAGGATTAAGTTTCAGGAGGTGATTCTGAGTCGATGGGTTTCTTC
ACGAGAGAGGAGTGACCAAGACGATCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC227271 representing NM_001145394
Red=Cloning site Green=Tags(s)


MAEAVFHAPKRKRVRVYETYESPLPIFGQDHGPLEKFIFRAEMINNVIVRNAEDIEQLYGKGYFGKGI
LSRSRPSFTISDPKLVAKWDMKTNMPIITSKRYQHSVEWAAELMRRQGQDESTVRRILKDYTKPLEHPP
VKRNEEAQVHDKLNSGMVSNMEGTAGGERPSVRGLHHEDGSQHIGLLHPGDRGPDHEYVLVEEAECAMSE
REAAPNEELVQRNRLICRRNPYRIFEYLQSLLEAFLVYALGCLSIYYEKEPLTIVKLWKAFTVVQPTF
RTTYMAYHYFRSGWVPKVGLKYGTDLLLYRKGPPFYHASYSVIELVDDHFEGLRRLPSWKSLAALS
RSVSVNYSKELMLCYLIKPMSTMDKEMESPECMKRIKVQEIVLSRWSSRRSDQDDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8057_g09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001145394

ORF Size: 1218 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145394.1](#), [NP_001138866.1](#)

RefSeq ORF: 1221 bp

Locus ID: 80746


UniProt ID: [Q8NCE0](#)

Cytogenetics: 3p25.2


MW: 46.9 kDa

Gene Summary: This gene encodes one of the subunits of the tRNA splicing endonuclease. This endonuclease catalyzes the first step in RNA splicing which is the removal of introns. Mutations in this gene have been associated with pontocerebellar hypoplasia type 2. A pseudogene has been identified on chromosome 4. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2009]

Product images:



Circular map for RC227271



Western blot validation of overexpression lysate (Cat# [LY428858]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227271 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).