

## Product datasheet for **RC227281L4V**

### Periostin (POSTN) (NM\_001135934) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Periostin (POSTN) (NM_001135934) Human Tagged ORF Clone Lentiviral Particle
Symbol:	POSTN
Synonyms:	OSF-2; OSF2; PDLPOSTN; PN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001135934
ORF Size:	2337 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227281).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001135934.1</a>
RefSeq ORF:	2340 bp
Locus ID:	10631
UniProt ID:	<a href="#">Q15063</a>
Cytogenetics:	13q13.3
Protein Families:	Druggable Genome, Secreted Protein
MW:	86.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a secreted extracellular matrix protein that functions in tissue development and regeneration, including wound healing, and ventricular remodeling following myocardial infarction. The encoded protein binds to integrins to support adhesion and migration of epithelial cells. This protein plays a role in cancer stem cell maintenance and metastasis. Mice lacking this gene exhibit cardiac valve disease, and skeletal and dental defects. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]