

Product datasheet for **RC227291L3V**

TCHP (NM_001143852) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TCHP (NM_001143852) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TCHP
Synonyms:	TpMs
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001143852
ORF Size:	1494 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227291).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001143852.1
RefSeq Size:	3111 bp
RefSeq ORF:	1497 bp
Locus ID:	84260
UniProt ID:	Q9BT92
Cytogenetics:	12q24.11
MW:	61.1 kDa



[View online »](#)

Gene Summary:

Tumor suppressor which has the ability to inhibit cell growth and be pro-apoptotic during cell stress. Inhibits cell growth in bladder and prostate cancer cells by a down-regulation of HSPB1 by inhibiting its phosphorylation. May act as a 'capping' or 'branching' protein for keratin filaments in the cell periphery. May regulate K8/K18 filament and desmosome organization mainly at the apical or peripheral regions of simple epithelial cells (PubMed:15731013, PubMed:18931701). Is a negative regulator of ciliogenesis (PubMed:25270598).[UniProtKB/Swiss-Prot Function]