

## Product datasheet for RC227332

### SPATA3 (NM\_139073) Human Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
Product Name: SPATA3 (NM\_139073) Human Tagged ORF Clone  
Tag: Myc-DDK  
Symbol: SPATA3  
Synonyms: TSARG1  
Mammalian Cell Selection: Neomycin  
Vector: pCMV6-Entry (PS100001)  
E. coli Selection: Kanamycin (25 ug/mL)  
ORF Nucleotide Sequence: >RC227332 representing NM\_139073  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAAGAAGGTCAAGAAGAAAAGGTCAAGAGGCCAGACGCCACCGAGACTCCACCTCCCAGCATGCTAGCT  
CCAATTCCACCTCTCAGCAGCCTAGCCCTGAATCCACACCACAGCAGCCTAGTCTGAATCCACACCACA  
GCAGCCTAGCCCTGAATCCACACCACAGCATTCCAGCCTTGAACCACCTCCCGGAGCCAGCATTCCAA  
GCCCTTCAGCACCCGAAATCCGCCGCTCCTCTTGCTGCCTTTATCTCCAGATGCTAACGTGAAGGCAG  
CCCCTCAATCCAGGAAAGCAGGGCCTCTGATTCGCGCCGCCCGCATTCTGCTCCTGTGCCACTTGCCC  
CTGCAGCTCCGCTTGCTGGCGTCTGCTGGGGCTATGCCATAGCCGCATCTTCGATGTCTCTGCCTCGG  
GACTGGCAGATGGCGCCAGGGAGAGGACTCCCCAACCTGCTCACCTTCTACAGAAAATCTTCAAGAAAAC  
CCTCCAGTCATCGTAACGCGTGTCTCCAAGCCCTCGGAACTGTGGCTGTGGCTCTGGGGCTCTAGGAG  
CTGCCTACTACATCAC

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227332 representing NM\_139073  
Red=Cloning site Green=Tags(s)

MKKVKKRSEARRHRDSTSQHASSNSTSQQPSPSTPQQPSPSTPQQPSPSTPQHSSLETTSRQPAFQ  
ALPAPEIRRSSCCLLSPDANVKAAPQSRKAGPLIRAGPHSCSCATCPCSSACWRRLGLCHSRIFDVLPR  
DWQMAPGRGLPNLLTFYRKSSRKPSHRNACPPSPRNCGCGSGGSRSCLLHH


**SGP**TRRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_139073

**ORF Size:** 576 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_139073.3](#), [NP\\_620712.2](#)

**RefSeq ORF:** 579 bp


**Locus ID:** 130560

**UniProt ID:** [Q8NHX4](#)

**Cytogenetics:** 2q37.1

MW: 20.7 kDa

Product images:



Circular map for RC227332