

Product datasheet for **RC227334**

BAP31 (BCAP31) (NM_001139457) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BAP31 (BCAP31) (NM_001139457) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BAP31
Synonyms:	6C6-AG; BAP31; CDM; DDCH; DXS1357E
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC227334 representing NM_001139457 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGCCGAGGCGTCCTCCTTGGTGCCCTGGCACTGCTCTCCCGAAGAACGCCTTTCAGTTAAAC
GGGCGTCGAAATCTCGGGCTTCTGGGCAGGGATCGTCGGGAGAGGCCGCTCTGGACGTGTTGACACA
CGTGCTGGAGGGGCAGGAAACAAGCTCACATCTCCTGTGGAAACCTTCTAGCAACAGGATGAGTCTG
CAGTGGACTGCAGTTGCCACCTCCTCTATGCGGAGGTCTTTGTTGTGTTGCTTCTGCATTCCCTTCA
TTTCTCCTAAAAGATGGCAGAAGATTTCAAGTCCCGGCTGGTGGAGTTGTTAGTGTCTATGGCAACAC
CTTCTTTGTGGTTCTCATTGTCACTCTGTGCTTTGGTCATCGATGCCGTGCGCGAAATTCGGAAGTAT
GATGATGTGACGAAAAGGTGAACCTCCAGAACAATCCCGGGGCCATGGAGCACTTCCACATGAAGCTTT
TCCGTGCCAGAGGAATCTACATTGCTGGCTTTTCTTGTGCTGTCTTCTGCTTAGACGCCTGGT
GACTCTCATTTTCGAGCAGGCCACGCTGCTGGCCTCCAATGAAGCCTTTAAAAGCAGGCCGAGAGTGCT
AGTGAGGCGGCCAAGAAGTACATGGAGGAGAATGACCAGCTCAAGAAGGGAGCTGCTGTTGACGGAGGCA
AGTTGGATGTCGGGAATGCTGAGGTGAAGTTGGAGGAAGAGAACAGGAGCCTGAAGGCTGACCTGCAGAA
GCTAAAGACGAGCTGGCCAGCACTAAGCAAAAAGCTAGAGAAAGCTGAAAACAGGTTCTGGCCATGCCG
AAGCAGTCTGAGGGCCTACCAAGGAGTACGACCGCTTGTGGAGGAGCACGCAAAGCTGCAGGCTGCAG
TAGATGGTCCCATGGACAAGAAGGAAGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC227334 representing NM_001139457
Red=Cloning site Green=Tags(s)


MGAEASSSWCPGTALPEERLSVKRASEISGFLGQGSSGEAALDVLTHVLEGAGNKLTSSCGKPSSNRMSL
 QWTAVATFLYAEVFFVLLLCIPFISPKRWQKIFKSRLVELLVSYGNTFFVVLIVLVLVIDAVREIRKY
 DDVTEKYNLQNNPGAMEHFHMKLFRAQRNLYIAGFSLLSFLLRRLVTLISQQATLLASNEAFKKQAESA
 SEAAKKYMEENDQLKKGAAVDGGKLDVGNAEVKLEENRSLKADLQKLKDELASTKQKLEKAENQVLAMR
 KQSEGLTKEYDRLLLEHAKLQAAVDGPMDKKEE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg8004_d03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001139457

ORF Size: 939 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001139457.2](#), [NP_001132929.1](#)

RefSeq ORF: 942 bp

Locus ID: 10134

UniProt ID: [P51572](#)


Cytogenetics: Xq28

Protein Families: Druggable Genome, Transmembrane


MW: 34.6 kDa

Gene Summary: This gene encodes a member of the B-cell receptor associated protein 31 superfamily. The encoded protein is a multi-pass transmembrane protein of the endoplasmic reticulum that is involved in the anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi and in caspase 8-mediated apoptosis. Microdeletions in this gene are associated with contiguous ABCD1/DXS1375E deletion syndrome (CADD5), a neonatal disorder. Alternative splicing of this gene results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 16. [provided by RefSeq, Jan 2012]


Product images:



Circular map for RC227334



Western blot validation of overexpression lysate (Cat# [LY427962]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227334 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified BCAP31 protein (Cat# [TP327334]). The protein was produced from HEK293T cells transfected with BCAP31 cDNA clone (Cat# RC227334) using MegaTran 2.0 (Cat# [TT210002]).