

Product datasheet for **RC227357A1V**

Human NCR3 (NM_001145466) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human NCR3 (NM_001145466) AAV Particle
Tag: Myc-DDK
Symbol: NCR3
Synonyms: 1C7; CD337; LY117; MALS; NKp30
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC227357 representing NM_001145466
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCTGGATGCTGTTGCTCATCTTGATCATGGTCCATCCAGGATCCTGTGCTCTCTGGGTGTCCCAGC
CCCCTGAGATTCGTACCCTGGAAGGATCCTCTGCCTTCCTGCCCTGCTCCTTCAATGCCAGCCAAGGGAG
ACTGGCCATTGGCTCCGTACGTGGTCCGAGATGAGGTGGTTCAGGGAAGGAGGTGAGGAATGGAACC
CCAGAGTTCAGGGGCCGCTGGCCCCACTTGCTTCTCCCGTTTCCTCCATGACCACCAAGGCTGAGCTGC
ACATCCGGGACGTGCGAGGCCATGACGCCAGCATCTACGTGTGCAGAGTGGAGGTGCTGGGCCTTGCTGT
CGGGACAGGGAATGGGACTCGGCTGGTGGTGGAGAAAGAACATCCTCAGCTAGGGGCTGGTACAGTCTCTC
CTCCTTCGGGCTGGATTCTATGCTGTACGCTTTCTCTCTGTGGCCGTGGGCAGCACCGTCTATTACCAGG
GCAATATGCCAATCTACTCTCTCCGATTCCCCCACTC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227357 representing NM_001145466
Red=Cloning site Green=Tags(s)

MAWMLLLILIMVHPGSCALWVSQPPEIRTLLEGSSAFLPCSFNASQGRLAIGSVTWFRDEVVPGKEVRNGT
PEFRGRPLASSRFLHDHQAELHIRDVRGHDASIYVCRVEVLGLGVGTGNGTRLVVEKEHPQLGAGTVL
LLRAGFYAVSFLSVAVGSTVYYYQGYAKSTLSGFPQL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_001145466
ORF Size:	531 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001145466.1</u>
RefSeq ORF:	534 bp
Locus ID:	259197
UniProt ID:	<u>O14931</u>
Cytogenetics:	6p21.33
MW:	19.1 kDa