

Product datasheet for **RC227360**

ZNF577 (NM_001135590) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ZNF577
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC227360 representing NM_001135590
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGAAAAATGCCACGATTGTAATGTCTGTGAGGAGAGAGCAAGGCAGTTCTTCAGGGGAGGGGTCATTGT
CATTGAAGATGTGGCTGTGGCTTACCAGGGAGGAGTGGCAGTTTTTGGACCACTCTCAGAAGGTCTT
GTACAAGGAAGTAATGTTGGAGAACTACATCAACCTAGTATCAATAGGGTATCGAGGCACCAAGCCAGAT
TCGCTCTTCAAGTTGGAGCAAGGAGAACCCCGAGGATAGCAGAAGGAGCAGCCACAGTCAAACTGTGC
CAGGAGGGAAACCATGAATGCAGTGTGTGCGGGAGAGCCTTCTCCAGGAAAGCACAGCTTATTCACCA
TCAGAGAACTGAAAGAGGAGAGAAAACCCCATGGATGTGGTGAATGTGGGAAAACATTCATGAGGAAGATT
CAGCTCACTGAGCATCAGAGAACTCACACAGGAGAGAAAGCCCATGAATGTAGTGAATGTGGAAAAGCCT
TCTCCAGAAAGTCACAGCTCATGGTCCATCAGAGAACCCATACAGGAGAGAAAACCTACAGATGCAGCAA
ATGCGGAAAAGCCTTCAGCCGGAAGTGCCGGCTCAATAGACATCAGCGATCACATACTGGAGAGAACTC
TATGGGTGCAGTGTGTGGGAAAGCCTTTCTCAGAAGGCATACCTCACTGCACACAGAGACTTCACA
CAGGAGATAAGCCTTATAAATGCAGTGATTGTGGAAGAACCTTCTATTTTAAGTCAGACCTGACCAGACA
TCAGAGGATTCATACGGGAGAGAAAACCTTATGAATGTAGTGAGTGTGAAAAAGCCTTTAGAAGCAAGTCG
AAGCTCATTGAGCATCAGCGTACTCACACTGGAGAGAGACCATATTCATGCAGAGAATGTGGCAAGCCT
TTGCCACATGTGAGTCTCATTAAACATGAGAAAACCTCACATAAGAGAGACAGCCATAAATTCAGTGAC
GGTGGAGAAAACCTTCCCTCAAGGAGTCATACCTCCTTATACATGAGCGAACTCATACAAGAGCAAAAGACT
GTGAACACAGTACCTATAGAAATGCCTTCCCTCAGGAACCCCGCCATTGTTAAACAAGAGTGAGCGCCTAG
TGGGCAGAAATGTAGTATTGTGGAACAACCTTTTCCAAGAAATCAAGCCTTTGTAGTTAATCAGGAATT
TGAACAGAGAATAAGCCTCACAAATGAAGTGAATGTGGCCCCATCAGTAATAAATTATCTTGTATCTT
ACAGATATTGTATCAGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC227360 representing NM_001135590
Red=Cloning site Green=Tags(s)


MKNATIVMSVRREQSSGEGSLSFEDVAVGFTREEWQFLDQSQKVL YKEVMLENYINLVSIGYRGTKPD
SLFKLEQGEPPGIAEGAAHSQICPGGKPHCEVCGRAF SRKAQLIQHQRTERGEKPHGCGECKGTFMRKI
QLTEHQRTHTGEKPHECCECGKAFSRKSQLMVHQRTHTGEKPYRCSKCGKAFSRKCRNLNRHQRSHTGEKL
YGCVSCGKAFSQKAYLTAHQRLHTGDKPYKSCDCGRTFYFKSDLTRHQRJHTGEKPYECSECEKAFRSKS
KLIQHQRTHTERPYPSCRECGKAFAHMSVL IKHEKTHIRETAINSLTVEKPPSSRSHTSLYMSEL IQEQKT
VNTVPIEMPSSSGTPPLLNKSERLVGRNVVIVEQPFPFNQAFVFNQEFEQRISLTNEVNVAPSVINYILYL
DTPVSE

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8059_a10.zip

Restriction Sites: SqaI-MluI

Cloning Scheme:



ACCN: NM 001135590

ORF Size: 1278 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001135590.1](#), [NP_001129062.1](#)

RefSeq ORF: 1281 bp

Locus ID: 84765

UniProt ID: [Q9BSK1](#)


Cytogenetics: 19q13.41

Protein Families: Transcription Factors

MW: 48.2 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC227360