

Product datasheet for **RC227364**

UNC93A (NM_001143947) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UNC93A (NM_001143947) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UNC93A
Synonyms:	dj366N23.1; dj366N23.2; Unc-93A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC227364 representing NM_001143947 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACAGAAGTCTAAGGAACGTCCTTGTGGTTTCCTTTGGGTTCTGCTTCTTTACAGCCTATGGAG
 GTCTGCAGAGCCTGCAGAGCAGCCTGTACAGCGAGGAGGCGCTGGGTGTACAGCGCTCAGCACCTCTA
 TGGAGGCATGCTCCTGTCCTCCATGTTCTCCACCGCTCCTCATCGAGAGGCTGGGCTGCAAGGGGACC
 ATCATCTCTCCATGTGTGGCTACGTGGCCTTCTCCGTGGCAACTTCTTCGCCAGCTGGTACACTTTGA
 TCCCCACCTCCATACTGCTGGGACTCGGGGCCGCCCGCTGTGGTCTGCACAGTGCACATACCTCACGAT
 CACGGGAAACACACATGCAGAGAAGCGGGAAAGCGTGGCAAAGACATGGTGAACCAATTTTGGCATC
 TTCTTCTCATATTCCAGTCATCCGGTGTGTGGGGCAACTTGATCTCATCGCTGGTATTTGGCCAGACTC
 CCAGCCAAGGGAGTGGTGTCTGGCTGTCTGATGATAGCTGCGTTCTCCAACCCATACGAGATGTTCA
 GCGGGAAGTGAAGGAGAGAAGAAATCAGTACCTTTCTGGTCCACTTTACTGTGCACTTTCAAGCTATAT
 AGAGATAAACGCTGTGCCTCTTAATTCTGTGCCGTGTACAGTGGATTGCAGCAAGGATTCCTCTCCA
 GCGAATACACAAGGTCTATGTCACCTGCACCTGGGCATCCAGTTCGTGGCTACGTGATGATCTGCTT
 CTCGGCCACTGACGCGCTGTGCTCCGTGTTGTATGAAAGGTCTCGCAGTACACGGGCAGGGCTGTGCTG
 TACGTGCTGGGCGCGGTGACCCAGTGTCTGCATGATTGCCCTACTGCTGTGGAGACCTCGTGCTGACC
 ATCTGGCAGTGTTCTTCTGATTCTCTGGCCTGTGGGCGTGGCAGATGCCGTCTGGCAGACACAAAACAA
 TGCTCTCTACGGCGTTCTGTTTGAGAAGAGCAAGGAAGCTGCCTTCGCCAATTACCGCCTGTGGAGGCC
 CTGGGCTTCGTATTGCTTCGGGTACAGCATGTTTTGTGCGTGCACGTCAAGCTCTACATTCTGCTGG
 GGGTCTGAGCCTGACCATGGTGGCGTATGGGCTTGTGGAGTGCCTGGAGTCCAAGAACCCGATCAGACC
 CCACGCTCCAGGACAGGTCAACCAGGCAGAGGATGAAGAAATACAAACAAAAATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC227364 representing NM_001143947
 Red=Cloning site Green=Tags(s)


MDRSLRNVLVVSFGFLLLTAYGGLQSLQSSLYSEEGLVLTALSTLYGGMLLSSMFLPPLLIERLGCKGT
 IILSMCGYVAFSVGNFFASWYTLIPTSILLGLGAAPLWSAQCTYLITGNTHAEKAGKRGKDMVNQYFGI
 FFLIFQSSGVWGNLISSLVFGQTPSQSGSVLAVLMIAAFLQPIRDVQRESEGEKKSVPFWSTLLSTFKLY
 RDKRLCLLILLPLYSLQQGFLSSEYTRSYVTCTLGIQFVGVMICFSATDALCSVLGKVSQYTGRAVL
 YVLGAVTHVSCMIALLLWRPRADHLAVFFVFSGLWGVADAVWQTQNNALYGVLFEKSKEAAFANYRLWEA
 LGFVIAFGYSMFLCVHVKLYILLGVLSLTMVAYGLVECVESKNPIRPHAPGVNQAEDEEIQTKM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001143947

ORF Size: 1245 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001143947.2](#)

RefSeq ORF: 1248 bp

Locus ID: 54346


UniProt ID: [Q86WB7](#)

Cytogenetics: 6q27

Protein Families: Transmembrane

MW: 45.5 kDa

Product images:



Circular map for RC227364