

## Product datasheet for **RC227373L4V**

### KPNA7 (NM\_001145715) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	KPNA7 (NM_001145715) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KPNA7
Synonyms:	IPOA8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001145715
ORF Size:	1548 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227373).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145715.1</a>
RefSeq ORF:	1551 bp
Locus ID:	402569
UniProt ID:	<a href="#">A9QM74</a>
Cytogenetics:	7q22.1
MW:	56.8 kDa


[View online »](#)

**Gene Summary:**

The transport of molecules between the nucleus and the cytoplasm in eukaryotic cells is mediated by the nuclear pore complex (NPC), which consists of 60-100 proteins. Small molecules (up to 70 kD) can pass through the nuclear pore by nonselective diffusion while larger molecules are transported by an active process. The protein encoded by this gene belongs to the importin alpha family, and is involved in nuclear protein import, but exhibits different nuclear localization signal binding specificity compared to other members of the family. A pseudogene of this gene has been defined on chromosome 5. [provided by RefSeq, Jul 2016]