

Product datasheet for **RC227385L2V**

ANKRD30B (NM_001145029) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ANKRD30B
Synonyms:	NY-BR-1.1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001145029
ORF Size:	4176 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC227385).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001145029.1](#), [NP_001138501.1](#)

RefSeq ORF: 4179 bp

Locus ID: 374860

Cytogenetics: 18p11.21

MW: 157.9 kDa

