

Product datasheet for **RC227395A1V**

Human NCR3 (NM_001145467) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: NCR3
Synonyms: 1C7; CD337; LY117; MALS; NKp30
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC227395 representing NM_001145467
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCTGGATGCTGTTGCTCATCTTGATCATGGTCCATCCAGGATCCTGTGCTCTCTGGGTGTCCAGC
CCCCTGAGATTCGTACCCTGGAAGGATCCTCGCCTTCCTGCCCTGCTCCTCAATGCCAGCCAAGGGAG
ACTGGCCATTGGCTCCGTACGTGGTTCGAGATGAGGTGGTTCAGGGAAGGAGGTGAGGAATGGAACC
CCAGAGTTCAGGGCCGCTGGCCCCACTTGCTTCTTCCCGTTCTCCATGACCACCAGGCTGAGCTGC
ACATCCGGGACGTGCGAGGCCATGACGCCAGCATCTACGTGTGCAGAGTGGAGGTGCTGGGCCTTGGTGT
CGGGACAGGGAATGGGACTCGGCTGGTGGTGAGAAAGAACATCCTCAGCTAGGGGCTGGTACAGTCCTC
CTCCTTCGGGCTGGATTCTATGCTGTGAGCTTTCTCTGTGGCCGTGGGACGACCGTCTATTACCAGG
GCAAATGCCACTGTCACATGGGAACACACTGCCACTCCTCAGATGGGCCCGAGGAGTGATTCCAGAGCC
CAGATGTCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227395 representing NM_001145467
Red=Cloning site Green=Tags(s)

MAWMLLLILIMVHPGSCALWVSQPPEIRTLEGSSAFLPCSFNASQGRILAIGSVTWFRDEVVPGKEVRNGT
PEFRGRPLASSRFLHDHQAELHIRDVRGHDASIYVCRVEVLGLGVGTGNGTRLVVEKEHPQLGAGTVL
LLRAGFYAVSFLSAVGSTVYYQKCHCHMGTHCHSSDGPGRVPIEPRCP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_001145467
ORF Size:	570 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001145467.1
RefSeq ORF:	573 bp
Locus ID:	259197
UniProt ID:	O14931
Cytogenetics:	6p21.33
MW:	20.5 kDa