

Product datasheet for **RC227430**

NSP3 (SH2D3C) (NM_001142533) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NSP3 (SH2D3C) (NM_001142533) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NSP3
Synonyms:	CHAT; NSP3; PRO34088; SHEP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC227430 representing NM_001142533
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACCGAGCGGTGCAGCCTGTGGAGTGCCTGTGCGCCGCCCTGCTGCTTCTACCGCGCTCCTTCG
 TGCAGGTGCAGTTCTCCAAGGAGAAGTACATCCTGGACTCATCGCCAGAGAACTCCACAAGGAATTGGA
 GGAGGAGCTCAAACCTCAGCAGCACGGATCTCCGACCCATGCCTGGTACCATGGCCGCATCCCCGAGAG
 GTCTCGGAGACCTTGGTACAACGCAACGGCGACTTCTCATCCGGGACTCACTACCAGCCTGGCGGACT
 ATGTGCTCACGTGCCGCTGGCGCAACAGGCCCTTGCACTTCAAGATCAACAAGGTGGTGGTGAAGGCAGG
 CGAGAGCTACACACACATCCAGTACCTGTTTGGCAGGAGAGCTTTGACCAGTGCCCGCCCTCGTGCGC
 TATCATGTGGGCAGCCGCAAGGCTGTGTCAGAGCAGAGTGGTCCATCATCTACTGCCCGGTGAACCGCA
 CCTTCCCACTGCGCTACCTCGAGGCCAGCTATGGCCTGGGACAGGGGAGTAGCAAGCCTGCTAGCCCCGT
 CAGCCCCCTCAGGCCCAAGGGCAGCCACATGAAGCGGGCAGCGTCACCATGACCATGGGCTCACTGCT
 GACAAGGTCAACCGCAGCGATGGCTGCCCCACCACTACGTCGCTGCCCGCCCTCGGGACTCCATCCGCA
 GCTGTGCCCTCAGCATGGACCAGATCCCAGACCTGCACTCACCCATGTCGCCCATCTCCGAGAGCCCTAG
 CTCCCCTGCCTACAGCACTGTAACCCGTGTCCATGCCGCCCTGCAGCCCTTCTGCCACAGCATTGGCT
 GCCTCCCCTGTGCGCCGCGTTCAGTGAAGCCAGCTGTGTCCCGGAAGTGCCCAAAGACCCATGGGG
 AGTCAGACAAGGGCCCCACACCAGCCCCCTCCACACCCCTTGGCAAGGCCTCCCCGTACCATCACTCAG
 CAGCTACAGTGACCCGGACTCTGGCCACTACTGCCAGCTCCAGCCTCCCGTGGTGGCAGCCGAGAGTGG
 GCAGCGACTGAGACCTCCAGCCAGCAGGCCAGGAGCTATGGGGAGAGGCTAAAGGAAGTGCAGAAAATG
 GGGCCCCGAAGGGGACTGGGGCAAGACCTTACAGTCCCCATCGTGGAAAGTCACTTCTTCTTCAACCC
 GGCCACCTTCCAGTCACTACTGATCCCCAGGGATAACCGGCCACTGGAGGTGGGCCCTTCTGCGCAAGGTC
 AAGGAGCTGCTGGCAGAAGTGGATGCCCGGACGCTGGCCCGCATGTACCAAGGTGGACTGCCTGGTTG
 CTAGGATACTGGGCGTTACCAAGGAGATGCAGACCCTAATGGGAGTCCGCTGGGGCATGGAAGTGCAC
 CCTCCCCATGGCCGGCAGCTACGCCTAGACCTGCTGGAAGGTTCCACACCATGTCCATCATGCTGGCC
 GTGGACATCCTGGGCTGCACCGGCTCTGCGGAGGAGCGGGCAGCGCTGCTGCACAAGACCATTAGCTGG
 CGGCCGAGCTGCGGGGACTATGGGCAACATGTTAGCTTCGCGGCGGTATGGGTGCCCTGGACATGGC
 CCAGATTTCTCGGCTGGAGCAGACATGGGTGACCTGCGGCAGCGACACACAGAGGGTGCATCCTGTAC
 GAGAAGAAGCTCAAGCCTTTTCTCAAGAGCCTCAACGAGGGCAAAGAAGGCCCGCCGCTGAGCAACACCA
 CGTTTCTCATGTGCTGCCCTCATCACCTGTGGAGTGTGACTCGGCCCCACCAGAGGGCCCTGAGCC
 CTGGGGCAGCACGGAGCACGGCTGGAGGTGGTGTGGCTCACCTGGAGGCCGCCGCACAGTGGCACAC
 CACGGAGGCCGTACCACACCAATGCTGAAGTCAAGCTGCAGGGGTTCCAGGCCCGCCGGAGCTCCTGG
 AGGTGTTACAGCACGGAGTTCAGATGCGCCTTCTCTGGGGCAGTCAAGGTGCCAGCAGCAGCCAGGCCCG
 GCGCTATGAGAAGTTCGACAAGGTCCTCACTGCCCTGTCCACAAGTGGAACTGCTGTCCGCTCCAGC
 GAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227430 representing NM_001142533
 Red=Cloning site Green=Tags(s)

```


MTERCSLWSALSAAACCFYRGSFVQVQFSKEKYILDSSPEKLHKELEELKLSSTDLRSHAWYHGRIPRE
VSETLVQRNGDFLIRDLSLTDYVLTTCRWRNQALHFKINKVVVKAGESYTHIQYLFEQESFDHVPALVR
YHVGSRKAVSEQSGAIIYCPVNRFTPLRYLEASYGLGQGSSKSPASPVSPSPGKSHMKRRSVTMTDGLTA
DKVTRSDGCPSTSLRPRRDSIRSCALSMQIPDLHSPMSPISESPSPAYSTVTRVHAAPAAPSATALP
ASPVARRSSESQLCPGSAPKTHGESDKGPHTSPSHTLGKASPSPLSSYSDPDSGHYCQLQPPVRSREW
AATETSSQARSYGERLKELENGAPEGDWGTKFTVPIVEVTSSFNATFQSLLIIPRDNRPLEVGLLRKV
KELLAEDVARTLARHVTKVDCLVARILGVTKEMQTLMGVRWGMELLTLPHGRQLRDLLERFHTMSIMLA
VDILGCTGSAEERAALLHKTQLAAELRGTMGNMFSFAAVMGALDMAQISRLEQTWVTLRQRHTEGAILY
EKKLPFLKSLNEGKEGPPLSNTTFPHVLPITLLECDSPPEGPEPWGSTEHGVEVLAHLEAARTVAH
HGGLYHTNAEVKLQGFQARPELLEVFSTEFQMRLLWGSQGASSSQARRYEKFDKVLTALESHKLEPAVRSS
EL
  
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001142533

ORF Size: 2106 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142533.1](#), [NP_001136005.1](#)

RefSeq ORF: 2109 bp

Locus ID: 10044

UniProt ID: [Q8N5H7](#)


Cytogenetics: 9q34.11

Protein Families: Druggable Genome

MW: 77.2 kDa

Gene Summary: This gene encodes an adaptor protein and member of a cytoplasmic protein family involved in cell migration. The encoded protein contains a putative Src homology 2 (SH2) domain and guanine nucleotide exchange factor-like domain which allows this signaling protein to form a complex with scaffolding protein Crk-associated substrate. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]

Product images:



Circular map for RC227430