

## Product datasheet for **RC227446L3V**

### **GALP (NM\_001145546) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GALP (NM_001145546) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GALP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145546
ORF Size:	147 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227446).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145546.1</a> , <a href="#">NP_001139018.1</a>
RefSeq ORF:	150 bp
Locus ID:	85569
UniProt ID:	<a href="#">Q9UBC7</a>
Cytogenetics:	19q13.43
Protein Families:	Secreted Protein
MW:	5.1 kDa

[View online »](#)

**Gene Summary:**

This gene encodes a member of the galanin family of neuropeptides. The encoded protein binds galanin receptors 1, 2 and 3 with the highest affinity for galanin receptor 3 and has been implicated in biological processes involving the central nervous system including hypothalamic regulation of metabolism and reproduction. A peptide encoded by a splice variant of this gene, termed alarin, has vasoactive properties, displays antimicrobial activity against *E. coli*, and may serve as a marker for neuroblastic tumors.[provided by RefSeq, Nov 2014]