

Product datasheet for **RC227449A1V**

Human Methionine Sulfoxide Reductase A (MSRA) (NM_001135670) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: Methionine Sulfoxide Reductase A
Synonyms: PMSR
Mammalian Cell: None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC227449 representing NM_001135670
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTCTCGGCCACCCGAGGGCTTGCCAGCTCCTCCTCCACAGCCTTTTCCCGTCCCGAGGATGG
GCAACTCGGCCTCGAACATCGTCAGCCCCAGGAGGCCTTGCCGGGCCGGAAGGAACAGACCCCTGTAGC
GGCCAAACATCATGTCAATGGCAACAGAACGTGCAACCTTTCCAGAGGGAACACAGATGGCTGTATTT
GAAAAAAGTGGCCATGCAGAAGTCGTCGAGTGGTGTACCAGCCAGAACACATGAGTTTGGAGAACTGC
TCAAGTCTTCTGGGAGAATCACGACCCGACCAAGGTATGCGCCAGGGGAACGACCATGGCACTCAGTA
CCGCTCGGCATCTACCCGACCTCTGCCAAGCAAATGGAGGCAGCCCTGAGCTCCAAGAGAACTACCAA
AAGGTTCTTTTCAGAGCACGGCTTCGGCCCCATCACTACCGACATCCGGGAGGGACAGACTTTCTACTATG
CGGAAGACTACCACCAGCAGTACCTGAGCAAGAACCCCAATGGCTACTGCGGCCTTGGGGGCACCGCGT
GTCCTGCCAGTGGGTATTAATAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC227449 representing NM_001135670
Red=Cloning site Green=Tags(s)

MLSATRRACQLLLLHSLFPVPRMGNSASNIQVSPQALPGRKEQTPVAAKHVNGNRTVEPFPEGTMAMV
EKTGHAEEVVRVYQPEHMSFEELLKVFWENHDPQGMROGNDHGTQYRSAYPTSAKQMEALSSKENYQ
KYLSEHGFPIITDIREGQTFYYAEDYHQYLSKNPNGYCGLGGTGVSCPVGIIK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_001135670
ORF Size:	585 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001135670.1
RefSeq ORF:	588 bp
Locus ID:	4482
UniProt ID:	Q9UJ68
Cytogenetics:	8p23.1
MW:	21.6 kDa