

Product datasheet for **RC227505**

STRA6 (NM_001142619) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STRA6 (NM_001142619) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	STRA6
Synonyms:	MCOPCB8; MCOPS9; PP14296
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC227505 representing NM_001142619
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCGTCCCAGCCAGCAGGGAACCCAGACCTCCCCGGGGCCACAGAGGACTACTCCTATGGCAGCTGGT
ACATCGATGAGCCCCAGGGGGCGAGGAGCTCCAGCCAGAGGGGAAGTGCCCTCCTGCCACACCAGCAT
ACCACCCGGCCTGTACCACGCCTGCCTGGCCTCGCTGTCAATCCTTGTGTGCTGCTCCTGGCCATGCTG
GTGAGGGCGCCAGCTCTGGCCTGACTGTGTGCGTGGCAGGCCCGCCTGCCAGGCCCGGGCAGTGC
CTGCTGCTGTTTTCATGGTCTCCTGAGCTCCCTGTGTTTGTGCTCCCCGACGAGGACGCATTGCCCTT
CCTGACTCTCGCCTCAGCACCCAGCAAGATGGGAAAAGTGGGCTCCAAGAGGGGCCTGGAAGATACTG
GGACTGTTCTATTATGCTGCCCTCTACTACCCTCTGGCTGCCTGTGCCACGGCTGGCCACACAGCTGCAC
ACCTGCTCGGCAGCACGCTGTCTGGGCCACCTTGGGGTCCAGGTCTGGCAGAGGGCAGAGTGTCCCCA
GGTGCCCAAGATCTACAAGTACTACTCCCTGCTGGCCTCCCTGCCTCTCCTGCTGGGCCTCGGATTCCTG
AGCCTTTGGTACCTGTGACAGCTGGTGAGAAGCTTCAGCCGTAGGACAGGAGCAGGCTCCAAGGGGCTGC
AGAGCAGTACTCTGAGGAATATCTGAGGAACCTCCTTTCAGGAAGAAGCTGGGAAGCAGCTACCACAC
CTCCAAGCATGGCTTCCTGTCTGGGCCCGCTGTGCTTGAGACTGCTACTACTCCACAGCCAGGA
TTCCATCTCCCGTGAAGCTGGTGTCTTTCAGCTACTGACAGGGACGGCCATTTACCAGGTGGCCCTGC
TGCTGCTGGTGGGCGTGGTACCCACTATCCAGAAGGTGAGGGCAGGGGTACCACGGATGTCTCCTACCT
GCTGGCCGGCTTTGGAATCGTGTCTCCGAGGACAAGCAGGAGGTGGTGGAGCTGGTGAAGCACCATCTG
TGGGCTGGAAGTGTGTACATCTCAGCCTTGGTCTTGTCTGCTTACTCACCTTCTGGTCTGATGC
GCTCACTGGTGACACACAGGACCAACCTTCGAGCTCTGCACCGAGGAGCTGCCCTGGACTTGAGTCCCTT
GCATCGGAGTCCCCATCCCTCCCGCCAAGCCATATTCTGTTGGATGAGCTTCAGTGCCTACCAGACAGCC
TTTATCTGCCTTGGGCTCCTGGTGCAGCAGATCATCTTCTTCTGGGAACCACGGCCCTGGCCTTCTGG
TGCTCATGCCTGTGCTCCATGGCAGGAACCTCCTGCTTTCGGTTCCTGGAGTCCCTGTGGCCCTTCTG
GCTGACTTTGGCCCTGGCTGTGATCCTGCAGAACATGGCAGCCATTGGGTCTTCTGGAGACTCATGAT
GGACACCCACAGCTGACCAACCGCGAGTGTCTATGCAGCCACCTTCTTCTTCCCCCTCAATGTGC
TGGTGGGTGCCATGGTGGCCACCTGGCGAGTGTCTCTGCCCCTTACAACGCCATCCACCTTGGCCA
GATGGACCTCAGCCTGTGCCACCGAGAGCCGCACTCTCGACCCCGGCTACTACAGTACCGAACTTC
TTGAAGATTGAAGTCAGCCAGTCGCATCCAGCCATGACAGCCTTCTGCTCCCTGCTCCTGCAAGCGCAGA
GCCTCCTACCCAGGACCATGGCAGCCCCCAGGACAGCCTCAGACCAGGGGAGGAAGACGAAGGGATGCA
GCTGCTACAGACAAAGGACTCCATGGCCAAGGGAGCTAGGCCCGGGGCCAGCCGGGCAGGGCTCGCTGG
GGTCTGGCCTACACGCTGCTGCACAACCAACCTGCAGGTCTTCCGCAAGACGGCCCTGTTGGGTGCCA
ATGGTGCCAGCCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227505 representing NM_001142619
 Red=Cloning site Green=Tags(s)


MSSQPAGNQTSPGATEDYSYGSWYIDEPQGGEEELQPEGEVPSCHTSIPPGLYHACLASLSILVLLLLAML
 VRRRQLWPDCVRGRPLPRPRAVPAAVFMVLLSSLCLLLPDEDALPFLTLASAPSDGKTEAPRGAWKIL
 GLFYAALYYPLAACATAGHTAAHLLGSTLSWAHLGVQVWQRAECPQVPKIYKYISLLASLPLLLGLGFL
 SLWYPVQLVRSFRRRTGAGSKGLQSSYSEEYLRNLLCRKKLGSSYHTSKHGFLSWARVCLRHCIYTPQPG
 FHLPLKLVLSATLTGTAIYQVALLLVGVVPTIQKVRAGVTTDVSYLLAGFGIVLSEDKQEVVELVKHHL
 WALEVCYISALVLSCLLTFVLVLRSLVTHRTNLRALHRGAALDLSPLHRSPHPSRQAIFCWMSFSAYQTA
 FICLGLLVQQIIFFLGTALAFVLMPVLHGRNLLFRSLESSWPFWLTLALAVILQNMAAHWVFLETHD
 GHPQLTNRRLYAATFLLFPLNVLVGAMVATWRVLLSALYNAIHLQMDLSLLPPRAATLDPGYYTYRNF
 LKIEVSQSHPAMTAFCSLLLQAQSLLPRTMAAPQDSLPGEEDEGMQLLQTKDSMAKGARPGASRGRARW
 GLAYTLLHNPTLQVFRKTALLGANGAQP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8071_g07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142619

ORF Size: 1974 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142619.2](#)

RefSeq ORF: 1977 bp

Locus ID: 64220

UniProt ID: [Q9BX79](#)


Cytogenetics: 15q24.1

Protein Families: Transmembrane

MW: 72.4 kDa

Gene Summary: The protein encoded by this gene is a membrane protein involved in the metabolism of retinol. The encoded protein acts as a receptor for retinol/retinol binding protein complexes. This protein removes the retinol from the complex and transports it across the cell membrane. Defects in this gene are a cause of syndromic microphthalmia type 9 (MCOPS9). Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Dec 2008]

Product images:



Circular map for RC227505