

Product datasheet for **RC227512A1V**

Human HNRNPCL2 (NM_001136561) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: HNRNPCL2

Synonyms: HNRNPCP5

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC227512 representing NM_001136561
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCAGCAACGTTACCAACAAGATGGATCCTCACTCCGTGAACCTCCCGTGTGTTTCATTGGGAATCTCA
ACACTCTTGTGTCAAGAAATCTGATGTGGAGGCGATCTTTTCCAAGTATGGCAAAATTCGGGGCTGCTC
TGTTTCATAAGGGCTTTGCCTTCGTTCAATATGATAAGGAGAAAAATGCCCGGGCTGCTGTAGCAGGAGAG
GATGGCAGAATGATTGCTAGCCAGGTTGCAGTTATTAACCTGGCTGCAGAGCCAAAAGTGAACCGAGGAA
ACGCAGGTGTGAAACGATCCGCAGCGGAGATGTACGGCTCCTCTTTGACTTGGACTATGGCTTTCAACG
GGATTATTATGATGGGATGTACAGTTTCCAGCACGTGTACCTCCTCCTCCTCCATTGCTCTGGCTGTA
GTGCCCTCGAAACGCCAGCGCATATCAGGAAACACCTCACGAAGGGGCAAAAGTGGCTTCAATTCTAAGA
GTGAAAGCGGGGATCTTCCAAGTCTGAAAGCTGAAAGGAGATGACCTTCAGGCCATTAAAGCAGGAGTT
GACCCAGATAAAACAGAAAGTGGATTCTCTCTGGAACCTGGAAAAATGAAAAGGAACAGAGCAAA
CAAGAGGTAGAGGTGAAAAATGCTAAGTCAGAAGAGGAGCAGAGCAGTAGCTCCATGAAGAAAGATGAGA
CTCATGTAAAGATGGAGTCTGAGGGGGTGCAGAAGACTCTGCTGAGGAGGGGGACCACTGGATGATGA
TGATAATGAAGATCAGGGGGACAACAGCTTCATTTGATCAAGAATAATGAAAAGATGCTGAGGAAGGA
GAGGATAACAGAGACAGCAACCAATGGCCAGGATGACTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC227512 representing NM_001136561 Red=Cloning site Green=Tags(s) MASVNTNKMDPHSVNSRVFIGNLNTLVVKKSDVEAIFSKYGKIAGCSVHKGFQVQYDKEKNARAAVAGE DGRMIASQVAVINLAAEPKVNRRGNAGVKRSAAEMYGSSFDLDYGFQRDYYDGMYSFPAVPPPPPIALAV VPSKRQRISGNTSRRGKSGFNSKSGKRGSSKSGKLKGDDLQAQKQELTQIKQKVDLLENLEKIEKEQSK QEVEVKNKSEEEQSSSSMKKDETHVKMESEGAEDSAEEDPLDDDDNEDQGDNLHLIKNNEKDAEEG EDNRDSTNGQDDSD TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001136561
ORF Size:	879 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001136561.1</u>
RefSeq ORF:	882 bp
Locus ID:	440563
UniProt ID:	<u>B2RXH8</u>
Cytogenetics:	1p36.21
MW:	31.9 kDa