

## Product datasheet for **RC227556**

### Choline Acetyltransferase (CHAT) (NM\_001142934) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Choline Acetyltransferase (CHAT) (NM_001142934) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CHAT
Synonyms:	CHOACTASE; CMS1A; CMS1A2; CMS6
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC227556 representing NM\_001142934  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCAGCAAAACTCCCAGCAGTGAGGAGTCTGGGCTGCCAAACTGCCCGTGCCCCCGCTGCAGCAGA  
 CCCTGGCCACGTACCTGCAGTGCATGCGACACTTGGTGTCTGAGGAGCAGTTCAGGAAGAGCCAGGCCAT  
 TGTGCAGCAGTTTGGGGCCCTGGTGGCCTCGGCGAGACCCTGCAGCAGAAACTCCTGGAGCGGCAGGAG  
 AAGACAGCCAACCTGGGTGTCTGAGTACTGGCTGAATGACATGTATCTCAACAACCGCCTGGCCCTGCCTG  
 TCAACTCCAGCCCTGCCGTGATCTTTGCTCGGCAGCACTTCCCTGGCACCAGTACCAGCTGAGGTTTGC  
 AGCCAGCCTCATCTCTGGTGTACTCAGCTACAAGGCCCTGCTGGACAGCCACTCCATTCCCAGTACTGT  
 GCCAAAGGCCAGCTGTCAGGGCAGCCCTTTGCATGAAGCAATACTATGGGCTCTTCTCCTCTACCGGC  
 TCCCCGGCCATACCCAGGACACGCTGGTGGCTCAGAACAGCAGCATCATGCCGGAGCCTGAGCACGTCAT  
 CGTAGCCTGCTCAATCAGTTCTTTGTCTTGGATGTTGTCATTAATTTCCGCCGTCTCAGTGAGGGGGAT  
 CTGTTCACTCAGTTGAGAAAGATAGTCAAAATGGCTTCCAACGAGGACGAGCGTTTGCCTCCAATTGGCC  
 TGCTGACGTCTGACGGGAGGAGCGAGTGGGCCGAGGCCAGGACGGTCTCGTGAAGACTCCACCAACCG  
 GGACTCGCTGGACATGATTGAGCGCTGCATCTGCCTTGTATGCCTGGACGCGCCAGGAGCGTGGAGCTC  
 AGCGACACCCACAGGGCACTCCAGCTCCTTACGGCGGAGGCTACAGCAAGAACGGGGCCAATCGCTGGT  
 ACGACAAGTCCCTGCAGTTTGTGGTGGGCCGAGACGGCACCTGCGGTGTGGTGTGCAAACTCCCCATT  
 CGATGGCATCGTCTGGTGCAGTGCAGTGCATCTGCTCAAGCACATGACGCAGAGCAGCAGGAAGCTG  
 ATCCGAGCAGACTCCGTCAGCGAGCTCCCCGCCCGGGAGGCTGCGGTGAAATGCTCCCCGAAATTC  
 AAGGCCACTTAGCCTCCTCGGCAGAAAACCTTCAACGAATAGTAAAGAACCTTGACTTCATTGTCTATAA  
 GTTTGACAACTATGGGAAAACATTCATTAAGAAGCAGAAATGCAGCCCTGATGCCTTCATCCAGGTGGCC  
 CTCCAGCTGGCCTTCTACAGGCTCCACCGAAGACTGGTGCCACCTACGAGAGCGCTCCATCCGCCGAT  
 TCCAGGAGGGACGCGTGGACAACATCAGATCGGCCACTCCAGAGGCACTGGCTTTTGTGAGAGCCGTGAC  
 TGACCACAAGGCTGCTGTGCCAGCTTCTGAGAAGCTTCTGCTCCTGAAGGATGCCATCCGTGCCAGACT  
 GCATACACAGTCATGGCCATAACAGGGATGGCCATTGACAACCACCTGCTGGCACTGCGGGAGCTGGCC  
 GGGCCATGTGCAAGGAGCTGCCGAGATGTTTCATGGATGAAACCTACCTGATGAGCAACCGGTTTGTCT  
 CTCCTACTGCCAGGTGCCACAACCACGGAGATGTTCTGCTGCTATGGTCTGTGGTCCCAATGGGTAT  
 GGTGCCCTGCTACAACCCCGCCAGAGACCATCTTTTCTGCATCTCTAGCTTTCACAGCTGCAAAGAGA  
 CTCTTTCTAGCAAGTTTGCAAAAGCTGTGGAAGAAAGCCTCATTGACATGAGAGACCTCTGCAGTCTGCT  
 GCCGCTACTGAGAGCAAGCCATTGGCAACAAAGGAAAAGCCACGAGGCCAGCCAGGGACACCAACT

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>RC227556 representing NM\_001142934  
 Red=Cloning site Green=Tags(s)

MAAKTPSSEESGLPKLPVPLQQTALATYLQCMRHLVSEEQFRKSQAIIVQFQFVAPGGLGETLQQLLERQE  
 KTANWVSEYWLNDMYLNNRLALPVNSSPAVIFARQHFPGTDDQLRFAASLISGLVSYKALLDSHSIPTDC  
 AKGQLSGQPLCMKQYYGLFSSYRLPGHTQDTLVAQNSSIMPEPEHVIVACCNQFFVLDVVINFRRLSEGD  
 LFTQLRKIVKMASNEDERLPIIGLLTSDGRSEWAEARTVLVKDSTNRDSLDMIERICLVCLDAPGGVEL  
 SDTHRALQLLHGGYSKNGANRWYDKSLQFVVGRTGTCGVVCEHSPFDGIVLVQTEHLLKHMTQSSRKL  
 IRADSVSELPAAPRRRLRWKCSPEIQGHLASSAEKLQRIKLNLDIVYKFDNYGKTFIKKQKSPDAFIQVA  
 LQLAFYRLHRRLLVPTYESASIRRFQEGRVDNIRSATPEALAFVRAVTDHKAAPASEKLLLLKDAIRAQT  
 AYTVMAITGMAIDNHLALRELARAMCKELPEMFMDETYLSNRFLVSTSQVPTTTEMFCCYGPVVPNGY  
 GACYNPQPETILFCISSFHCKETSSSKFAKAVEESLIDMRDLCSLLPPTESKPLATKEKATRPSQGHQP


**SGP**TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001142934

ORF Size: 1890 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001142934.1</a> , <a href="#">NP_001136406.1</a>
<b>RefSeq ORF:</b>	1893 bp
<b>Locus ID:</b>	1103
<b>UniProt ID:</b>	<a href="#">P28329</a>
<b>Cytogenetics:</b>	10q11.23
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Glycerophospholipid metabolism
<b>MW:</b>	70.2 kDa
<b>Gene Summary:</b>	This gene encodes an enzyme which catalyzes the biosynthesis of the neurotransmitter acetylcholine. This gene product is a characteristic feature of cholinergic neurons, and changes in these neurons may explain some of the symptoms of Alzheimer's disease. Polymorphisms in this gene have been associated with Alzheimer's disease and mild cognitive impairment. Mutations in this gene are associated with congenital myasthenic syndrome associated with episodic apnea. Multiple transcript variants encoding different isoforms have been found for this gene, and some of these variants have been shown to encode more than one isoform. [provided by RefSeq, May 2010]