

Product datasheet for RC227558

C15orf61 (NM_001143936) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: Myc-DDK
Symbol: C15orf61
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC227558 representing NM_001143936
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGCCCTGAGGAGGGCCACGAGGTCGCGCTCCGCCTGCTGCTGTAGGCCGTGGCCTCGCGG
 CCGCCGCCGCCCAAGCCAGCGCCTCGGAGGTGCTGACGCGGCATCTGCTGCAGCGCGCCTGCCGCA
 CTGGACCTCCTTCTGCGTGCCCTACAGCGCCGTCGCAACGACCAAGTTCGGCCTCTCGCACTTCAACTGG
 CCGGTGCAGGGCGCAACTACCACGTCTGCGCACCGGCTGCTTCCCTTCATCAAGTACCACTGCTCCA
 AGGCTCCCTGGCAGGACCTGGCCCGGCAGAACCGCTTCTTACGCGCTCAAGGTCGTCAACCTCGGTAT
 TCCAACTTTATTATAGGACTTGCTCCTGGTTATTTGCCAGAGTCACAGAGACTGTGCATACCAGTTAT
 GGACCCATAACAGTTTATTTCTCAATAAAGAAGATGAAGGTGCCATGTAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227558 representing NM_001143936
 Red=Cloning site Green=Tags(s)

MEALRRAHEVALRLLLCRPWASRAAAPKPSASEVLTRHLLQRRLPHWTSFCVPYSAVRNDQFGLSHFNW
 PVQGANYHVLRTGCFPIKYHCSKAPWQDLARQNRFFALKVNVNLIPTLLYGLGSWLFARVTVHTSY
 GPITVYFLNKEDEGAMY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8006_f11.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001143936

ORF Size: 471 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001143936.2](#)

RefSeq ORF: 474 bp


Locus ID: 145853

UniProt ID: A6NNL5
 Cytogenetics: 15q23
 Protein Families: Transmembrane
 MW: 17.9 kDa

Product images:



Circular map for RC227558



Western blot validation of overexpression lysate (Cat# [LY428408]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227558 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).