

Product datasheet for **RC227558A1V**

Human C15orf61 (NM_001143936) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: C15orf61
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC227558 representing NM_001143936
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGGCCCTGAGGAGGGCCACGAGGTCGCGCTCCGCTGCTGTGTAGGCCGTGGGCCTCGCGCG
CCGCCCGCCCAAGCCAGCGCCTCGGAGGTGCTGACGCGGCATCTGCTGAGCGGCGCCTGCCGCA
CTGGACCTCCTTCTGCGTGCCCTACAGCGCCGTCGCAACGACCAGTTCGGCCTCTCGCACTTCAACTGG
CCGGTGCAGGGCGCAACTACCACGTCCTGCGCACCGGCTGCTTCCCTTCATCAAGTACCACTGCTCCA
AGGCTCCCTGGCAGGACCTGGCCCGGAGAACCGCTTCTTACGGCGCTCAAGGTCGTCAACCTCGGTAT
TCCAACCTTTATTATATGGAAGTGGCTCCTGGTTATTTGCCAGAGTCACAGAGACTGTGCATACCAAGTTAT
GGACCCATAACAGTTTATTTTCTCAATAAAGAAGATGAAGGTGCCATGTAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227558 representing NM_001143936
Red=Cloning site Green=Tags(s)

MEALRRAHEVALRLLCRPWASRAAARPKPSASEVLTRHLLQRRLPHWTSFCVPYSAYVRNDQFGLSHFNW
PVQGANYHVLRTGCFPFIKYHCSKAPWQDLARQNRFFTALKVVNLGIPTLLYGLGSWLFARVETETVHTSY
GPITVYFLNKEDEGAMY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001143936



ORF Size:	471 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001143936.1
RefSeq ORF:	474 bp
Locus ID:	145853
UniProt ID:	A6NNL5
Cytogenetics:	15q23
MW:	17.9 kDa