

Product datasheet for **RC227618L3V**

HBS1L (NM_001145207) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HBS1L (NM_001145207) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HBS1L |
| Synonyms: | EF-1a; eRF3c; ERFS; HBS1; HSPC276 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001145207 |
| ORF Size: | 1896 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC227618). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001145207.1 |
| RefSeq ORF: | 1899 bp |
| Locus ID: | 10767 |
| UniProt ID: | Q9Y450 |
| Cytogenetics: | 6q23.3 |
| MW: | 69.9 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the GTP-binding elongation factor family. It is expressed in multiple tissues with the highest expression in heart and skeletal muscle. The intergenic region of this gene and the MYB gene has been identified to be a quantitative trait locus (QTL) controlling fetal hemoglobin level, and this region influences erythrocyte, platelet, and monocyte counts as well as erythrocyte volume and hemoglobin content. DNA polymorphisms at this region associate with fetal hemoglobin levels and pain crises in sickle cell disease. A single nucleotide polymorphism in exon 1 of this gene is significantly associated with severity in beta-thalassemia/Hemoglobin E. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, May 2009]