

Product datasheet for RC227624L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

BANF1 (NM_001143985) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: BANF1 (NM_001143985) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: BANF1

Synonyms: BAF; BCRP1; D14S1460; NGPS

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC227624).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





st The last codon before the Stop codon of the ORF.

ACCN: NM_001143985

ORF Size: 267 bp





BANF1 (NM_001143985) Human Tagged Lenti ORF Clone - RC227624L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001143985.1</u>, <u>NP 001137457.1</u>

 RefSeq Size:
 1122 bp

 RefSeq ORF:
 270 bp

 Locus ID:
 8815

 UniProt ID:
 075531

Cytogenetics: 11q13.1 **MW:** 10.1 kDa


Gene Summary: The protein encoded by this gene was first identified by its ability to protect retroviruses from

intramolecular integration and therefore promote intermolecular integration into the host cell genome. The protein forms a homodimer which localizes to both the nucleus and cytoplasm and is specifically associated with chromosomes during mitosis. This protein binds to double stranded DNA in a non-specific manner and also binds to LEM-domain containing proteins of the nuclear envelope. This protein is thought to facilitate nuclear reassembly by binding with both DNA and inner nuclear membrane proteins and thereby recruit chromatin to the nuclear periphery. Alternative splicing results in multiple transcript variants encoding


the same protein.[provided by RefSeq, Jan 2009]



Product images:



Circular map for RC227624L1



Double digestion of RC227624L1 using Sgfl and Mlul