

Product datasheet for **RC227640**

Consortin (CNST) (NM_001139459) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Consortin (CNST) (NM_001139459) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CNST
Synonyms:	C1orf71; PPP1R64
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC227640 representing NM_001139459
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATGACAGCGATACTCCTACATATTATCTGCAAAATAGAACCAAGATGGATGTCATCTGGTGACA
 CGGTGAAAGGAGTGTGACCTGTCTGCCTTCTGCATCAGATGAAAATGAAAATCAGCTTGACGGGGACGG
 GCATGAGCATCTGACCAGCAGTGACAGTGGGAAAGCCCAAGTGTCTGAGCAGGACAGTCTCAAT
 AATAATGAAAGCTGCACATTGAGCTGCGAGGTGGCTGCAGGTGAGAACTCGCAAAACACCCCTTTGTGAAG
 CCTCCAGAGATGAACAGGCCTTCTGGGAAAGGACAAAAAAATTCCTGGAAAAAGAGTCCAAGAAGCAA
 AAAAGGGACTGCTAAGAAGATACCACCAGGACTTTTTTCAGGAGATATTGCACCTTAATGCAAGAAAA
 GTACTAAGCGCAGTCACATATGCTGTTGATGATGAAGAAGCTGCTGAAGTAAATGCTAATGAGCAGCCAG
 AGGCGCCAAAGCTTGTCTGCAGTCTCTGTTTTCACTTATACGAGGTGAAGTTGAGCAGTTGGATTCAAG
 AGCACTTCCCCTTGCCTTCATCAGATAGCAGAATCCTATTTCCAGGAGGAGACTATGAGAAAGCAATG
 AAATTCATTGAGCTAGAAGCATTGTATCATGAGCAATTGCTCGCAAACTTTCTGCCATTCAAGAACAGT
 GGGAAACAAAATGGAAAAGTGTGCAACACATACAGTTACGGCTCTAAGGAATTCAGAAAAGGGATTTAA
 TGGTGAAGATTTGAACGGCTTACGAAAATTTGTGCAACACATCAAGATCCTCTTTTATCCAAACATAAG
 ATAGCAGCTGTGGAGAAGTCCCAGGAAAGGAAATGCTCCACGCAGTTACTAGTGTCTGAAGATCCAAAGG
 AAGGAGGAGCTACCACCAAAGAGTCAGAGAGTAAAAGTGTCTCGGCACAGAGTCAAGTAAAGAAAGCCA
 ACATACAGTGGAGCCCTGGGGAGCAGTCCCTGCTGTATCAGATGGACGTGCAACAGATTCCCAAGC
 CTTTCGGTAAGTGCAGGAAAGGACCACATGGAGGAGCTGCTGACGCGCTGAAGCCACGTTAGCGCTCC
 ACACCCAGTCTCCGAGACAGCAGGGAGCCCTGCTGGCCAGACTCTTCTGAGGATGCTTGTGAGGATGA
 CAGTCGCTTGCAGCTGGCTCAAACAGAGCCCTGCCAGGATGTGGCCAGAATAGAGGGCATTGCTGAAGAC
 CCTAAGGTGTTTTCTTCCAGCAAGTCAAAGACAGAGCCGTTGATTTACCAGGCTGTGACCGTATACCTC
 CTGCATTGATTTCTGAGGGTAAATATTCACAGGCTCAGAGGAAAGAACTCCGTTTGCCACTTCGGGATGC
 TTCTGAGGCGTTGCCACAGACCAACTTGAGAACAATGAATTAATGAGCTGCAGCAGCCTGATCTTACA
 GACAGTGTGAAATACCCACAGGCGCAGGCTGACTCAGACGGTTCTGAGAATGTGCTCTGTGAAATA
 ATCAAATATCTGACTTAGGCATACTGCTCCAGAGGTGTGTATGGCCCCAGAGGAAAAGGGAGATAAAGA
 CGACCAACTCAACAAAGAAACAGAAGACTATTTGAACAGCCTTTTGAAGGATGTTTAAAGATACTGAA
 GATTCCCTTTCTATGAAGATAACCAAGACGACTCCGATCTCCTTCAAGATCTCTCTCTGAAGAAG
 CATCCTATAGTCTCCAGGAGAATCTGCCTTCTGATGAGAGCTGTCTTTCTTGTGATGATCTTCCAAAAG
 GATAGAGATTGCAGAGGTAATCAGAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC227640 representing NM_001139459
 Red=Cloning site Green=Tags(s)

MDDSDTPTYLLQIEPQDGCHPGDSVERSVCPLPSASDENENQLDGDGHEHLTSSDSAMGKPQVSEQDSL
 NNESTLSCEVAAGENSQNTLCEASRDEQAF LGDKKIPGKRSPRSKKGTAKKIPPGLFSGDIAPLMQEK
 VLSAVTYAVDDEEAAEVNANEQPEAPKLVLSLIRGEVEQLDSRALPLCLHQIAESYFQEEDYEKAM
 KFIQLERLYHEQLLANLSAIQEQQWETKWKTVQPHVTALRNSEKGFNGEDFERLTKICATHQDPLL SKHK
 IAAVEKSQERKQSTQLLVSEDPKEGGATTKESKTCCLGTESSKESQHTVEPLGSSPCCHQMDVQTDSPS
 LSVTAGKDHMEELLCSAEATLALHTQSSETAGSPSPGDSSEDACEDDSRLQLAQTEACQDVARIEGIAED
 PKVFLSSKSKTEPLISPGCDRIPPALISEGKYSQAQRKELRLPLRDASEALPTDQLENNELNELQQPDLT
 DSDGKSPQAQADSDGSENVLCGNNQISDLGILLPEVCMAPEEKGDKDDQLNKETEDYLNLSLLEGCKLKDTE
 DSLSYEDNQDDSDLLQLDSPEEASYSLQENLPSDESCLSLDDLAKRIEIAEVNQR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001139459

ORF Size: 1848 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001139459.1 , NP_001132931.1
RefSeq ORF:	1851 bp
Locus ID:	163882
UniProt ID:	Q6PJW8
Cytogenetics:	1q44
Protein Families:	Transmembrane
MW:	67.4 kDa
Gene Summary:	Targeting of numerous transmembrane proteins to the cell surface is thought to depend on their recognition by cargo receptors that interact with the adaptor machinery for anterograde traffic at the distal end of the Golgi complex. Consortin (CNST) is an integral membrane protein that acts as a binding partner of connexins, the building blocks of gap junctions, and acts as a trans-Golgi network (TGN) receptor involved in connexin targeting to the plasma membrane and recycling from the cell surface (del Castillo et al., 2010 [PubMed 19864490]). [supplied by OMIM, Jun 2010]