

Product datasheet for **RC227648L3V**

CNG1 (CNGA1) (NM_001142564) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CNG1 (CNGA1) (NM_001142564) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CNGA1 |
| Synonyms: | CNCG; CNCG1; CNG-1; CNG1; RCNC1; RCNCa; RCNCalpha; RP49 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001142564 |
| ORF Size: | 2277 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC227648). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001142564.1 , NP_001136036.1 |
| RefSeq ORF: | 2061 bp |
| Locus ID: | 1259 |
| UniProt ID: | P29973 |
| Cytogenetics: | 4p12 |
| Protein Families: | Druggable Genome, Ion Channels: Cyclic nucleotide gated, Transmembrane |
| MW: | 86.7 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is involved in phototransduction. Along with another protein, the encoded protein forms a cGMP-gated cation channel in the plasma membrane, allowing depolarization of rod photoreceptors. This represents the last step in the phototransduction pathway. Defects in this gene are a cause of retinitis pigmentosa autosomal recessive (ARRP) disease. Multiple transcript variants have been found for this gene. [provided by RefSeq, Oct 2019]