

Product datasheet for **RC227669A1V**

Human C17orf107 (NM_001145536) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Human C17orf107 (NM_001145536) AAV Particle
 Tag: Myc-DDK
 Symbol: C17orf107
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >RC227669 representing NM_001145536
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAAGGGGACCCCCAGCTCCCTGGACACCCTGATGTGGATCTACCACTTCCACAGCTCCACCGAGGTGG
 CCCTCCAGCCCCCGTCTCTTCCCTGGAACCTCCGTGGCCGAGCCCATGAATATCTGGAGCAGAG
 GTTCAGAGAGCTGAAGTCCCTGGAGCCACCCGAACGAAGATGCAGGGGATGCTGCCTGCCCGAAGCCC
 ACCCTGGGGCTGGTGTGAGAGAAGCCACAGCCAGCCTCGTGAGCTTCGGCACCACTTGTAGAGATCT
 CAGCCCTGTGGCTGCAGCAGGAGGCGCGCGACTAGACGGCAGCGGGGCCAGCCCAGACGGCAGGGA
 TCCGGGTGCAGCGCTGAGCCGGGTAGCCCAAGCCGAGGGCAGGGGTTTCGGCAAGCTGGGGCTGCGGTG
 GGCGCCAGCGCCCGGCTCCTAGTCCAGGGAGCATGGCTATGCCTGTGTGGACGGGTCTGCAGGGGTCTG
 CCTATTCTGCGACAGTCGCAACAGCAGCTAGGCCTCGGAATCCCCGGAGAACCAGGTGAGCTCAGGACA
 CGGGGTGAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227669 representing NM_001145536
 Red=Cloning site Green=Tags(s)

MKGTPSSLDLTMWIYHFHSSTEVALQPPLSSLELSVAAAHEYLEQRFRELKSLEPPEPKMQGMLPAPKP
 TLGLVLREATASLVSGTTLLEISALWLQEARLDGSAGAPDGRDPGAALSRVAQAAGQGVQAGAAGV
 GASARLLVQGAWLCLCGRGLQGSASFLRQSQQLGLGIPGEPVSSGHGVS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
 Serotype: AAV-2



[View online »](#)

ACCN:	NM_001145536
ORF Size:	570 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001145536.1</u>
RefSeq ORF:	573 bp
Locus ID:	100130311
UniProt ID:	<u>Q6ZR85</u>
Cytogenetics:	17p13.2
MW:	19.8 kDa