

Protein Sequence:	>RC227684 representing NM_001145463 Red=Cloning site Green=Tags(s) MHPAGLAAAAAGTPRLPSKRRIIPVSQPGMADPHQLFDDTSSAQSRGYGAQRAPGGLSYPAASPTPHAAFL ADPVSNMAMAYGSSLAAQGKELVDKNIDRFIPITKLKYYFAVDTMVYGRKLGLLFFPYLHQDWEVQYQQD TPVAPRFDVNAPDLYIPAMAFITYVLVAGLALGTQDRFSPDLLGLQASSALAWLTLEVLAILLSLYLVTV NTDLTTIDLVAFLGYKYVGMIGGVLMLLFGKIGYYLVLGWCCVAIFVFMFLLPGAVAHACNPSTLGGR GGRITRSGRCG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001145463
ORF Size:	873 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001145463.1</u>
RefSeq ORF:	876 bp
Locus ID:	90522
UniProt ID:	<u>Q5BJH7</u>
Cytogenetics:	19q13.2
MW:	31.1 kDa