

## Product datasheet for **RC227686A1V**

### Human P311 (NREP) (NM\_001142478) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human P311 (NREP) (NM\_001142478) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** P311  
**Synonyms:** C5orf13; D4S114; P311; PRO1873; PTZ17; SEZ17  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC227686 representing NM\_001142478  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

**ATGGTTTATTACCCAGAACTCTTTGTCTGGGTCAAGAACCAATTTCCAAACAAGGACATGGAGGGAA**  
**GGCTTCCTAAGGGAAGACTTCCTGTCCCAAAGGAAGTGAACCGCAAGAAGACGATGAGACAAACGCTGC**  
**CTCCCTGACTCCACTGGCAGCAGTGAAGTCCGCTCCCAAGAATCAGTTACCTCCACTTTTT**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC227686 representing NM\_001142478  
**Red**=Cloning site **Green**=Tags(s)

MVYYPELFWVVSQEPFPNKDMEGRLPKGRLPVPKEVNRKKNDETNAASLTPLGSELRSRPRISYLHFF

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001142478  
**ORF Size:** 204 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

**Stability:** AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

**RefSeq:** [NM\\_001142478.1](#)

**RefSeq Size:** 2081 bp

**RefSeq ORF:** 207 bp

**Locus ID:** 9315

**UniProt ID:** [Q16612](#)

**Cytogenetics:** 5q22.1

**MW:** 7.9 kDa