

## Product datasheet for **RC227693**

### PRKRA (NM\_001139517) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRKRA (NM_001139517) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PRKRA
Synonyms:	DYT16; HSD14; PACT; RAX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC227693 representing NM_001139517 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCAAAGCACACCTTTTGTGGTTTTTGCAGTTTGGGAAGATGATAACAGCTAAGCCAGGAAAAACAC  
CGATTCAGGTATTACACGAATACGGCATGAAGACCAAGAACATCCAGTTTATGAATGTAAAGATCTGA  
TGTGCAAATACACGTGCCACTTTCACCTTCAGAGTAACCGTTGGTGACATAACCTGCACAGGTGAAGGT  
ACAAGTAAGAAGCTGGCGAAACATAGAGCTGCAGAGGCTGCCATAAACATTTTGAAGCCAATGCAAGTA  
TTTGCTTTGCAGTTCCTGACCCCTTAATGCCTGACCCCTCCAAGCAACCAAGAACCAGCTTAATCCTAT  
TGGTTCATTACAGGAATTGGCTATTCATCATGGCTGGAGACTTCCTGAATATACCCTTTCCAGGAGGGA  
GGACCTGCTCATAAGAGAGAATATACTACAATTTGCAGGCTAGAGTCATTTATGGAACTGGAAAGGGG  
CATCAAAAAGCAAGCCAAAAGGAATGCTGCTGAGAAATTTCTTGCCAAATTTAGTAATATTTCTCCAGA  
GAACCACATTTCTTAAACAAATGTAGTAGGACATTCCTTAGGATGTACTTGGCATTTCCTTGAGGAATTCT  
CCTGGTAAAAGATCACTTACTGAAAAGAAGCCTCCTTAGTATTCCAAATACAGATTACATCCAGCTGC  
TTAGTGAATTGCCAAGGAACAAGGTTTTAATAAACATATTTGGATATAGATGAACTGAGCGCCAATGG  
ACAATATCAATGTCTTGCTGAACTGTCCACCAGCCCCATCACAGTCTGTATGGCTCCGGTATCTCCTGT  
GGCAATGCACAAAGTGATGCAGCTCACAAATGCTTTGCAGTATTTAAAGATAATAGCAGAAAGAAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC227693 representing NM\_001139517  
Red=Cloning site Green=Tags(s)

MQSTPFCGFCSLGKMITAKPGKTPIQVLHEYGMKTKNIPVYECERSDVQIHVPTFTFRVTVDITCTGEG  
 TSKKLAKHRAAEAAINILKANASICFVPDPLMPDPSKQPKNQLNPIGSLQELAIHHGWRLPEYTLSEQE  
 GPAHKREYTTICRLESFMETGKGASKKQAKRNAAEKFLAKFSNISPENHISL TNVVGHS LGCTWHSLRNS  
 PGEKINLLKRSLLSIPNTDYIQLLSEIAKEQGFNITYLDIDEL SANGQYQCLAELSTSPITVCHGSGISC  
 GNAQSDAAHNALQYLKIIAERK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8050\\_g05.zip](https://cdn.origene.com/chromatograms/mk8050_g05.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001139517


**ORF Size:** 906 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001139517.1, NP_001132989.1</u>
<b>RefSeq ORF:</b>	909 bp
<b>Locus ID:</b>	8575
<b>UniProt ID:</b>	<u>O75569</u>
<b>Cytogenetics:</b>	2q31.2
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	32.9 kDa
<b>Gene Summary:</b>	This gene encodes a protein kinase activated by double-stranded RNA which mediates the effects of interferon in response to viral infection. Mutations in this gene have been associated with dystonia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2008]

### Product images:



Circular map for RC227693



Western blot validation of overexpression lysate (Cat# [LY427977]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC227693 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).