

Product datasheet for **RC227709L4V**

JMJD5 (KDM8) (NM_001145348) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	JMJD5 (KDM8) (NM_001145348) Human Tagged ORF Clone Lentiviral Particle
Symbol:	JMJD5
Synonyms:	JMJD5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001145348
ORF Size:	1362 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227709).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145348.1 , NP_001138820.1
RefSeq ORF:	1365 bp
Locus ID:	79831
UniProt ID:	Q8N371
Cytogenetics:	16p12.1
MW:	50.7 kDa
Gene Summary:	This gene likely encodes a histone lysine demethylase. Studies of a similar protein in mouse indicate a potential role for this protein as a tumor suppressor. Alternatively spliced transcript variants have been described.[provided by RefSeq, Feb 2009]


[View online »](#)