

## Product datasheet for **RC227724L4V**

### **KIAA0141 (DELE1) (NM\_001142603) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KIAA0141 (DELE1) (NM_001142603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DELE1
Synonyms:	DELE; DELE1(L); KIAA0141
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142603
ORF Size:	1545 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227724).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001142603.1</a>
RefSeq Size:	4195 bp
RefSeq ORF:	1548 bp
Locus ID:	9812
UniProt ID:	<a href="#">Q14154</a>
Cytogenetics:	5q31.3
MW:	55.9 kDa



[View online »](#)

**Gene Summary:**

Essential for the induction of death receptor-mediated apoptosis through the regulation of caspase activation.[UniProtKB/Swiss-Prot Function]