

Product datasheet for **RC227735L3V**

ADAMDEC1 (NM_001145271) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAMDEC1 (NM_001145271) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAMDEC1
Synonyms:	M12.219
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145271
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227735).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145271.1 , NP_001138743.1
RefSeq Size:	2438 bp
RefSeq ORF:	1176 bp
Locus ID:	27299
UniProt ID:	O15204
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome, Secreted Protein
MW:	52.8 kDa



[View online »](#)

Gene Summary:

This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells. [provided by RefSeq, Jul 2008]