

Product datasheet for **RC227765**

C8orf45 (MCMDC2) (NM_001136160) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C8orf45 (MCMDC2) (NM_001136160) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C8orf45
Synonyms:	C8orf45
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC227765 representing NM_001136160
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCAAATCTAAAAATGAAAGAGGCGGCCCTCATCTATCTTGACAGAAGTGAGGCCTCCAAAAGTTTA
TAGATGATTGCAAGTACTACAATGATTCAAAACAAAGCTATGCTGTCTATCGATTCAAAATTTAATAAA
TCCTCTGATGTTGTTGAATTAGATGCTGAGCTTGAAATCACATTTTACACCAACCTTTAAAAGCTGCT
GAAGTCTTTCAATCAGTCTGTTTTATTGCTGTTAAGACTCTCTCATTAAATTGGACAATTGCAGACTGAAA
CGCAAATTAATATAGTGTGAAATTAACACATTTACCTCCCTGCCAAGTTATGGTCTTGATCTTTGTGA
GTTTCCACTTGATTATACATCTCAGAGATTTTATATGATGCAAGGAATTGTGATTGCAATGACAATA
ACCAAGTATACACAAGGGCAAGATTTCTTTGTTGAGTGAAGCATGCCCTCTTCAAAGGATTTTCAGT
ATATAAGAGTGCATGTGCTGGTGTACAGAATCTGCAACGATAAGAAAATGACTTTTTGTGAATCTATG
TGCATCTTCACTCAAGAAGACAGAAAATTTAGAGTACTTGGTGATAACAAAATAGTTGAAATAATTGCC
ACAAAGGCACTTCGTGCTTTTCAAGGATATTCTAACAACAGCCATTTAGGTTTCAATCACTTACAATTT
TCTAAGAGATGAATCAGTGAATAAAATGAATATAGGAAATGAATATAAAATTATTGGAATTCACCTG
TGTA AAAACCTCACA AACTGCTGTCTGTATAGAAGCAATAGCATAACTTTTTGTAATTCAAAAGTTCT
TCAGGAATTAGTGACA AACTTCAGGTGTCTCCTCTCCTTAACTTCCAGCTCATGCTGGAAGTTTACAGCA
TACTTGCCAAATCTTTGCATCACA AATACCCCTCCTGGGACTTACAATTTGCTCAAGCTCTGTTTGT
GATGAGTCTAGTACAGACA AACTGACCGTAAACAAGGAAGTGAAGATTGCCTGGATATTTAATTATAACA
AGTGATACTCTACTCATAGACAGGCTTCTGAATTTTAGCATAAACCTTGCCCCGTGGTATACGTCATC
TAGTCTCTACTGAAATTTTCCCACTCTATCCAGGAATAAGTATGGAAGTGGAGCAGTTAGCATTACAGGC
TGGCAGTGCTTTGCTAGCTAAAGGTGGTATCTGCTTTATAGGAGACTTGGCTTACACAAAAAAGATAAA
CTTGAACAGCTTCAAACAGTTCTGGAGAGCAGAAGCATCACAGTGTACATCCCAGGAAAAGAGTTTGGG
AGGATATTGATCAACAGATGACTTTTCCAGTTCAGTGCAGTTTTTGGTCTTTTGTGATGTGGATTCATC
TTCAAGGAGAAAATGCACAGAAAATCAACACTCTAATTGGTCAGATGGATTGCAGTTTGATTCCAGCTAAC
CTTGTGGAGGCATTTGGTTTATTGATTA AACTGCAACGAGTCATCTCCTGCCACCCATTTCTCCTACTG
TGCAACACACTTTGAACAAAGCCATTAATCCTGAAGGGCTGTTTTATGCGGCTTCTAGACAGTTCACAAC
TGAAGATTTTAAAAGCTATTGGCTTTTGCAAAGAATTTGAATGTAGAATTCAGCTTGAAGCAGAAAAGA
ATGACCCATGGCTATTATCTAGCAAGTCGAGAATCAGAACAGGCTCTGTATGTGGATCAAAGCTGTCAG
CATCTGCATTA AAAATATCTCGTTTTCTATCTGAAGCCATGCACGACTGAACTTAAGGAACAAAGTGCT
TAAAGAAGATGTGCTGATTGCAGCCTTACTATTTGAAACATCCCTCACATTGAAATATGGTAATGAATTA
AATTAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227765 representing NM_001136160
 Red=Cloning site Green=Tags(s)


MSNLKMKEAALIYLDRSGGLQKFIDDCYYNDSKQSYAVYRFKILINPSDVVELDAELGNHILHQPLKAA
 EVFQSVCFIAVKTL SLIGQLQTEQTQINI VLKLT HLPPLPSYGLDLCEFPDYTSQRFYMMQGIVIAMTTI
 TKYTQGARFLCSDEACPLSKGFQYIRVHVPGATESATIRNDFLCNL CASSLQEDRKFRVLGDKQIVEIIA
 TKALRAFQGYSNQPPFRFQSLTIFLRDESVNKMNI GNEYKIIIGIPTCVKTSQTAVCIEANSITFCNSKVP
 SGISDNFRCLLSLTSSSCWKFTAILANIFASQITPPGTYNLLKLC LLSLVQTTDRNKELEDCLDILIIT
 SDTLLIDRLLNFSINLVPRGIRHLVSTEIFPTLSRNKYGTGAVSIQAGSALLAKGGICFIGDLASHKDK
 LEQLQTVLESRSITVYIPGKKFGEDIDQMTFPVQCSFWSFVDVSDSSRRAQKINTLIGQMDCSLIPAN
 LVEAFGLLINCNESSPCHPFLPTVQHTLNKAINPEGLFYAASRQFTTEDFEKLLAFAKNLNVEFSLEAER
 MTHGYLASRRIRTGSVCGSKLSASALKYL VFLSEAHARLNLRN KVLKEDVLIALLFETSLTKYGNEL
 NY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8070_f06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001136160

ORF Size: 1896 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001136160.1](#), [NP_001129632.1](#)

RefSeq ORF: 1899 bp


Locus ID: 157777

Cytogenetics: 8q13.1

MW: 70.5 kDa

Gene Summary: Plays an important role in meiotic recombination and associated DNA double-strand break repair.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC227765