

Product datasheet for **RC227789L4V**

EFTUD2 (NM_001142605) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EFTUD2 (NM_001142605) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EFTUD2
Synonyms:	MFDGA; MFDM; SNRNP116; Snrp116; Snu114; U5-116KD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142605
ORF Size:	2811 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227789).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001142605.1 , NP_001136077.1
RefSeq Size:	4058 bp
RefSeq ORF:	2814 bp
Locus ID:	9343
UniProt ID:	Q15029
Cytogenetics:	17q21.31
Protein Pathways:	Spliceosome
MW:	105.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a GTPase which is a component of the spliceosome complex which processes precursor mRNAs to produce mature mRNAs. Mutations in this gene are associated with mandibulofacial dysostosis with microcephaly. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]